



Mueller Thumb Stabilizer



The Mueller Thumb Stabilizer orthosis provides firm support to limit CMC and MCP joint motion while still allowing free IP motion. This one size fits most splint is easy-to-fit and easy-to-wear to treat CMC arthritis, deQuervain's, soft tissue or ligament injuries and for Gamekeeper's or Skier's Thumb.

Features

- Fully cushioned for comfortable for all day and/or all night wear.
- Two flexible stays provide maximum support to limit MP and CMC joint motion.
- Breathable fabric is protected with AEGIS Microbe Shield Controls to control odors, unsightly stains, and product deterioration
- Cost and time effective, the Mueller Thumb Stabilization splints minimizes stocking and maximizes time with your patients.

Order Information

Measure circumference at the wrist.

One size fits 5.5" to 10.5" (13.9cm - 26.6cm)

Universal design fits right or left hand

Product Number: P2015-U

Watch our **instructional video on the Mueller Thumb Stabilizer** on our website:
<http://www.3pointproducts.com/product-videos-of-splints/mueller-thumb-stabilizer-video>



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com