



Protexgloves Original



With light compression and sun protection, Protexgloves Original full finger glove provides support and soothing comfort for patients with Arthritis, Psoriasis, Raynaud's Syndrome, Scleroderma, Ehlers Danlos Syndrome or Neuropathy. The unique blend of materials safeguard delicate skin without sacrificing dexterity.

Features

- These form-fitting gloves feature a unique blend of Supplex nylon and Lycra making them lightweight enough to perform most fine motor tasks easily
- Durable enough for gardening or cleaning chores, these gloves have 50+ UVA and UVB sun protection knit right in
- Excellent skin protection for patients with Scleroderma, EDS or Raynaud's
- Recommended for day and/or nighttime wear to decrease swelling and inflammation for less stiffness and improved movement
- Sold one pair (two gloves) per package

Order Information

Measure length from tip of middle finger to distal wrist crease and width across the MP joints.

Size	Length	Width
Small	6 1/2"	2 1/2"
Medium	7"	3"
Large	7 1/2"	3 3/4"

Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com