



# Radial Hinged Ulnar Deviation Splint



The Radial Hinged Ulnar Deviation Splint aligns each finger individually to improve the alignment of strongly-deviated fingers. The radial hinge allows free MP flexion and extension for improved hand function.

## Features

- Adjustable finger straps permit each finger to be individually aligned for improved functional positioning and maximum control.
- The hinge allows fingers to flex and extend while acting as a counterforce to the unwanted ulnar deviation.
- The Kydex shell is lined in soft felt and can be heated for custom fit by a healthcare provider.
- Developed in the clinic to improve strength and function in patients with advanced Rheumatoid Arthritis.
- All materials are latex-free.

## Order Information

Lay hand flat and measure width across the hand from the 2nd to the 5th MP (large knuckles).

Size	Right	Left	Fits
Sm	P2003-R2	P2003-L2	2 3/4" to 3"
Med	P2003-R3	P2003-L3	3" to 3 3/8"
Lrg	P2003-R4	P2003-L4	3 3/8" to 3 5/8"



**Call 888-378-7763**

Mon-Fri 8:30am - 5:00pm EST for ordering details

**[www.3pointproducts.com](http://www.3pointproducts.com)**