



3pp[®] Step Down[™] Splint



The 3pp[®] Step Down[™] Splint is an easy, effective way to restore combined DIP and PIP flexion following injury following injury or surgery. Static progressive stretch can be easily applied to tolerance for optimal treatment compliance.

Features

- Adjustable, lightweight wire splint allows the wearer to control the amount of force applied.
- The case hardened steel wire frame will not bend or break from tension.
- Adjusts from zero degrees to near full flexion.
- Combine with a 3pp Final Flexion Wrap to restore final degrees of flexion.
- Pads are included to apply progressive force at the DIP as needed. Latex-free.
- Sizing and fitting by a health care provider is recommended.

Order Information

Measure length from web space to fingertip and width at proximal phalanx. Fits right or left hand. If in between sizes, choose the larger size.

	Right or Left	Web Space to Tip
Sm	P1202-2	1 1/2" - 2" (3.8 - 5cm)
Med	P1202-3	1 3/4" - 3 1/8" (4.4 - 8.3cm)
Lrg	P1202-4	3" - 3 1/2" (7.6 - 8.9cm)



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com