



ThumSaver™ MP



The hand-based ThumSaver™ MP is an excellent choice to protect thumbs with CMC arthritis or ligament injuries. Great for ball players to wear inside their gloves. The ThumSaver MP is also suitable for Gamekeeper’s thumb, basal joint arthritis, Skier’s thumb, sprains and strains.

Features

- Reduce pain when pinching or gripping small objects. This thumb splint realigns the MP joint (middle joint) and stabilizes the base of the thumb to relieve the grinding pain of arthritis.
- Wearing the ThumSaver MP, the wearer can garden, write or grasp a doorknob without pain or fear of dropping things.
- Ideal for Gamekeeper’s thumb, the ThumSaver MP is thin enough to be worn comfortably under gloves.
- Felt pads are included to provide cushioning where needed.

Order Information

Measure the length from the IP joint (knuckle at the end of the thumb) to the wrist. Measure around the IP joint.

	Right	Left	Length	Circumference
Sm	P2001-R2	P2001-L2	3" to 3 1/2"	2 1/2" to 2 3/4"
Med	P2001-R3	P2001-L3	3 1/2" to 4"	2 3/4" to 3 1/8"
Lrg	P2001-R4	P2001-L4	4" to 4 1/2"	3 1/8" to 3 1/4"
X-Lrg	P2001-R5	P2001-L5	4 1/2" to 5"	3 1/4" to 3 1/2"



Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com