
PPE - NONWOVEN DISPOSABLE ISOLATION GOWNS



Techstyles® nonwoven disposable isolation gowns help promote healthier environments by providing quality infection control barriers for staff and visitors.

Proper selection of protective apparel is essential to successfully limiting exposure and reducing potential cross-contamination as rapidly emerging health threats continue to challenge facilities nationwide.

Our isolation gowns are available in a wide variety of styles and performance rated materials from standard fluid-resistance to impervious and three AAMI protection levels.



Better care starts with safety and comfort.

VALUECARE® NONWOVEN DISPOSABLE ISOLATION GOWNS









4550-27X

Half Back Detail

ValueCare® Disposable Isolation Gowns

CONSTRUCTION

CAT NO	SEAM	NECK	WRIST	BACK	MATERIAL	SIZE	SWEEP	COLOR	PACKAGING	QTY/CASE
4600-340	Sewn	Tie	Elastic	Full	SBPP	Universal	55.0"	Yellow	10/BG, 5/CS	50 EA
4600-350	Sewn	Tie	Elastic	Full	SBPP	X-Large	59.0"	Yellow	10/BG, 5/CS	50 EA
4600-370	Sewn	Tie	Elastic	Full	SBPP	3X-Large	69.0"	Yellow	10/BG, 5/CS	50 EA
4550-27X	Heat	Tie	Thumb Loop	Open	SMS	X-Large	37.5"	Yellow	10/BG, 10/CS	100 EA
4620-340	Sewn	Tie	Elastic	Full	SMS & SBPP	Universal	55.0"	Yellow	10/BG, 5/CS	50 EA
4620-350	Sewn	Tie	Elastic	Full	SMS & SBPP	X-Large	59.0"	Yellow	10/BG, 5/CS	50 EA
4540-140	Sewn	Tie	Elastic	Half	PE/SBPP	Universal	27.5"	White	10/BG, 5/CS	50 EA
4540-150	Sewn	Tie	Elastic	Half	PE/SBPP	X-Large	29.5"	White	10/BG, 5/CS	50 EA
4540-145	Sewn	Tie	Knit	Half	PE/SBPP	Universal	27.5"	White	10/BG, 5/CS	50 EA
4540-155	Sewn	Tie	Knit	Half	PE/SBPP	X-Large	29.5"	White	10/BG, 5/CS	50 EA
4640-140	Sewn	Tie	Elastic	Full	PE/SBPP	Universal	55.0"	White	10/BG, 5/CS	50 EA
4640-150	Sewn	Tie	Elastic	Full	PE/SBPP	X-Large	59.0"	White	10/BG, 5/CS	50 EA
4640-145	Sewn	Tie	Knit	Full	PE/SBPP	Universal	55.0"	White	10/BG, 5/CS	50 EA
4640-155	Sewn	Tie	Knit	Full	PE/SBPP	X-Large	59.0"	White	10/BG, 5/CS	50 EA
4550-25X	Heat	Apron	Thumb Loop	Open	PE	X-Large	37.5"	Blue	18/BX, 5 BX/CS	90 EA













SAFECARE® NONWOVEN DISPOSABLE ISOLATION GOWNS







SafeCare AAMI Gowns (46969-087, 46969-088, 46969-089)



SafeCare Non-AAMI Gowns (46969-097X, 46969-099, 46969-076)

SafeCare® AAMI Rated Disposable Isolation Gowns

AAMI			CONS	TRUCTION							
LEVEL	CAT NO	SEAM	NECK	WRIST	BACK	MATERIAL	SIZE	SWEEP	COLOR	PACKAGING	QTY/CASE
1	46969-087	Sonic	Velcro	Elastic	Full	SMS	Universal	57"	Yellow	10/BG, 10/CS	100 EA
1	46969-187	Sonic	Velcro	Elastic	Full	SMS	X-Large	61"	Blue	10/BG, 10/CS	100 EA
2	46969-088	Heat	Velcro	Elastic	Full	SMS	Universal	57"	Yellow	10/BG, 10/CS	100 EA
3	46969-089	Heat	Velcro	Elastic	Full	SMS	Universal	57"	Yellow	10/BG, 10/CS	100 EA

SafeCare® Disposable Isolation Gowns

~~	NICT		LOT	401
1 11	NS ⁻	I KI	11 11	11 11
	110	1110	/ 🔾 🗆	101

CAT NO	SEAM	NECK	WRIST	BACK	MATERIAL	SIZE	SWEEP	COLOR	PACKAGING	QTY/CASE
46969-075	Sonic	Velcro	Elastic	Full	SMS	Universal	57"	Yellow	10/BG, 10/CS	100 EA
46969-076	Sonic	Velcro	Elastic	Full	SMS	Universal	57"	Yellow	10/BG, 10/CS	100 EA
46969-099	Sonic	Velcro	Elastic	Half	PE/SMS	Universal	28.5"	White	10/BG, 5/CS	50 EA
46969-100	Sonic	Velcro	Elastic	Full	PE/SMS	Universal	57"	White	10/BG, 5/CS	50 EA
46969-098	Heat	Apron	Thumb Loop	Open	PE	Universal	35.5"	Yellow	100 Bulk	100 EA
46969-097X	Heat	Apron	Thumb Loop	Open	PE	X-Large	37.5"	Blue	1/BG, 15/BX, 5/CS	75 EA
46969097X10	Heat	Apron	Thumb Loop	Open	PE	X-Large	37.5"	Blue	10/BG, 10/CS	100 EA
46969-097XB	Heat	Apron	Thumb Loop	Open	PE	X-Large	37.5"	Blue	1/BG. 100/CS	100 EA

SafeCare® Pediatric Disposable Isolation Gowns

CONSTRUCTION									QTY/
CAT NO	SEAM	NECK	WRIST	BACK	MATERIAL	SIZE	COLOR	PACKAGING	CASE
46969-PSM	Sewn	Velcro	Elastic	Full	SMS	Small 2-4 Years	Yellow	10/BG, 5/CS	50 EA
46969-PMD	Sewn	Velcro	Elastic	Full	SMS	Medium 6-8 Years	Yellow	10/BG, 5/CS	50 EA
46969-PLG	Sewn	Velcro	Elastic	Full	SMS	Large 10-12 Years	Yellow	10/BG, 5/CS	50 EA




















LET ENCOMPASS BE YOUR GUIDE TO BETTER PROTECTION

AAMI Classification Levels of Barrier Performance

AAMI classified gowns are required to meet strict testing requirements for the entire isolation gown, including the front, back, sleeves and seams as these are all considered critical zones as shown.

Critical Zones



AAMI Level	Applicable Test	Required Result	Exposure Risk
1	AATCC 42 - Impact Penetration Test	<= 4.5 g	Minimal
2	AATCC 42 - Impact Penetration Test AATCC 127 - Hydrostatic Pressure Test	<= 1.0g >= 20 cm	Low
3	AATCC 42 - Impact Penetration Test AATCC 127 - Hydrostatic Pressure Test	<= 1.0g >= 50 cm	Moderate

Encompass Non-AAMI Fluid Resistance Classifications

Encompass offers a wide variety of gowns that are not AAMI classified due to construction characteristics. To assist you in properly identifying the fluid resistance capabilities of these items, please use the following category descriptions and look for the category icon next to the catalog numbers.

SYMBOL GLOSSARY

ICON	CLASSIFICATION	FLUID EXPOSURE RATING	SUGGESTED ENVIRONMENT
	Fluid Resistant	Designed for low exposure. The level of fluid resistance is minimal.	Protects against the possible occurrence of minimal splashing and spraying.
	Fluid Guard	Designed for moderate exposure. The level of fluid resistance is moderate.	Protects against the likely occurrence of moderate splashing and spraying.
	Fluid Barrier	Designed for high exposure. The level of fluid resistance is high.	Protects against the probable occurrence of spraying and limited soaking.
	Impervious	Designed for the highest exposure. Passes ASTM F1670, Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Synthetic Blood and ASTM F1671, Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.	Protects against soaking or strike through.

^{*}American Society of Testing and Materials (ASTM)

Final selection must consider the degree of exposure risk for various procedures and be determined by the purchaser. Encompass provides no warranty against transmission of blood-borne diseases. For more information about our protection products, please contact your Encompass sales representative.

