

# So quiet, patient and partner can rest well at night



Effective seal reduces undesired airflow and accompanying noise

The air vent is located on the frame, not the cushion, improving the way air is dispersed Unique vent hole design disperses air, gently and quietly, away from patient and partner

- Quiet—to meet the needs of patients and their partners.
- Quiet—thanks to valuable feedback from you, our therapy experts.



### Uniquely effective nasal pillows technology

Mirage Swift II performs well on the three essential criteria for effectiveness: fitting, seal and comfort.

The uniqueness of ResMed's nasal pillows technology is in the microscopic detail—you can't see it, but the user can feel the difference.

#### Adjustable angle-fitting just right

 Rotating pillow sleeve assembly enables the user to select the most suitable angle for a comfortable fit and effective seal.

#### Shape-every microscopic detail counts

- The shape of the Mirage Swift II pillows is the result of extensive anthropometric studies of a wide range of facial and nasal structures, from diverse racial backgrounds.
- Every tiny element of the shape and design helps optimize seal, fitting and comfort for a wide range of users.

#### Flexibility-it matters

- The nasal pillows, the sleeve on which they are located, and the "neck" connecting the two, are all supple and flexible.
- This enables the Mirage Swift II to adapt effectively to the user's facial characteristics and accommodates nasal or facial movement.

Unique shape and design optimize fitting, seal and comfort

Angle of pillows at the nostrils can be adjusted, enhancing fitting, seal and comfort


Flexible pillows and neck on a flexible sleeve helps system adapt and conform to the user





## Designed for real people with real needs

- First-time users, older patients, younger patients, those who are less dexterous ... Mirage Swift II is designed for real people in the real world.
- Fitting the Mirage Swift II is quick and intuitive. After first-time adjustments, if any, the Mirage Swift II can be fitted without clips or Velcro®.
- Easy to assemble, disassemble and clean, due to few, easy-to-attach-ordetach parts, as well as visible guides and indents.





Each system comes packed with three pillow sizes—one of them just right for each patient







## Nasal pillows technology leaps forward to enhance your business or service

- · Quick fitting
- · Reduced inventory (three pillow sizes packed with product)
- · More choice for your patients
- · Comprehensive patient support package



#### PRODUCT CODE

Europe	60513
UK	60516
France	60518
ANZ	60514
Asia Pacific	60515

The Mirage Swift II comes complete with a storage and travel bag, fitting guide, detailed user instructions, and a DVD containing fitting and care information for the patient

ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Brazil, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). Protected by patents: AU 724360, US 6561191. Other patents and design registrations pending. Mirage and Mirage Swift are trademarks of ResMed Ltd. Velcro is a trademark of Velcro Industries B.V. Specifications may change without notice. ©2006 ResMed Ltd. 1011018/3 06 09