



Mirage Swift™ II

NASAL PILLOWS SYSTEM

Quietly improving mask acceptance

Mirage Swift II is an effective, light and quiet alternative to nasal masks.

The system covers less of the face, which can facilitate therapy acceptance and long-term compliance. The pillows sit gently at the entrance of the nostrils, reducing contact with the face.

Mirage Swift II is designed to provide quiet therapy, so that patients and their partners can sleep well at night. It is also designed to minimize dead space and carbon dioxide rebreathing, which are important for effective therapy.

Mirage Swift II—from ResMed, the world leaders in mask technology.

Secure headgear design

enhances seal and stability, without the need for overtightening

Nasal pillows

sit flexibly and comfortably at the entrance of the nostrils, for secure seal during therapy

Soft and flexible pillow sleeve

helps the nasal pillows stay in place, while applying minimal pressure to the nose and face

Tube with rotating elbow and swivel

offers more choice of sleeping positions

ResMed logo

makes it easy to position the system the correct way

Legible size mark

is easy to read for reordering

Rotating pillow sleeve assembly

enables adjustment to choose the most comfortable angle of the nasal pillows toward the nostrils



ACTUAL SIZE



TECHNOLOGY HIGHLIGHT

Unique air vent for quiet therapy

- The air vent is located on the frame, not the cushion, improving the way air is dispersed
- Unique vent hole design disperses air, gently and quietly, away from patient and partner

PATIENT BENEFITS

- Quiet therapy—for patients and their partners
- Comfortable fit and effective seal—rotating pillow sleeve assembly enables users to position the nasal pillows at the angle most suitable for them
- Adapts effectively to the user's face and accommodates nasal or facial movement—the nasal pillows, the sleeve on which they are located, and the 'neck' connecting the two, are all supple and flexible
- Quick and intuitive fitting at home—after first-time adjustments (if any), the Mirage Swift II can be fitted without clips or Velcro™
- Freedom to read or watch TV with mask and glasses on at bedtime—light and compact system
- Choice of sleeping positions—the tube can be attached on either side, facing upward or downward

PROFESSIONAL BENEFITS

- Quick initial fitting—designed so that the system can be placed on the patient's head easily and fitted quickly, with minimal adjustments
- Designed to fit a wide range of facial and nasal structures, from diverse racial backgrounds—shape based on extensive anthropometric studies
- Reduced inventory—each system comes packed with three pillow sizes, with one to suit each patient
- More choice for your patients—an effective alternative to nasal masks
- Comprehensive patient support package—the Mirage Swift II comes complete with its own bag, fitting guide, detailed user instructions and a user-friendly DVD containing fitting and care information for the patient
- Can be recommended for first-time users and older, younger or less dexterous patients—ergonomic design with few parts and visible guides; easy to assemble, disassemble and clean

PRODUCT CODES

Mirage Swift II Nasal Pillows System (includes three pillow sizes: S, M, L)

Europe	60513
UK	60516
France	60518
ANZ	60514
Asia Pacific	60515

