

# RESMED

## Mask performance you can trust

Mirage Micro Nasal Mask provides personalised fit and comfort to more users (like you!) around the world.

Your mask comes with a comprehensive user support package, which includes travel bag and user guide with fitting and care instructions.



**Mirage Micro is the latest generation nasal mask from ResMed —leaders in mask technology.**

## Mirage Micro

NASAL MASK

### PROVEN SEAL, PERFECT FIT

Personalised fit and comfort

Quiet, sleek and stylish

Intuitively easy



Visit [www.resmed.com](http://www.resmed.com) for more information on Mirage Micro, including step-by-step fitting instructions.


ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Brazil, China, Finland, France, Germany, Hong Kong, India, Japan, Malaysia, Netherlands, New Zealand, Norway, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). Protected by patents: AU 710733, AU 741003, AU 766623, AU 775051, AU 777033, AU 785376, CA 2261790, EP 0956069, EP 1187647, EP 1187648, EP 1187649, EP 1187650, EP 1314445, EP 1479406, HK 1057714, JP 3686609, NZ 513052, NZ 526165, NZ 526166, NZ 526167, NZ 526168, US 6112746, US 6119693, US 6357441, US 6374826, US 6439230, US 6463931, US 6532961, US 6557556, US 6581594, US 6581602, US 6634358, US 6691708, US 6823865, US 6860269, US 6871649, US 6997188, US 7011090, US 7159587, US 7178527, US 7188620, US 7216647, US 7234466, US 7243651. Other patents pending. Protected by design registrations: AU 304309, CN 200530117106.1, EP 366844, JP 1264474, NZ 406357, US D532512, US D547439, US D545960. Others pending. Mirage and Mirage Micro are trademarks of ResMed Ltd and Mirage is registered in U.S. Patent and Trademark Office. Specifications may change without notice. © 2007 ResMed Ltd. 1012279/1 07 12

# RESMED

## Mirage Micro™

NASAL MASK

User Information



Proven seal,  
perfect fit



## Personalised fit and comfort

- **A fine-tuned fit to suit you**—MicroFit™ dial allows precise adjustments of the forehead support, setting the mask at the most comfortable position for you.
- **Proven seal**—Mirage™ dual-wall cushion technology reduces pressure on the nasal bridge and seals effectively.
- **A winning combination**—as you adjust the dial to your profile, the cushion works to seal effectively ... so there is no need to overtighten the headgear.
- **A range of sizes**—with five cushions and two frames to choose from, you'll find the perfect fit for you.

## Quiet, sleek and stylish

- **Quiet performance**—an array of small, built-in vents diffuses air gently and quietly.
- **Clear field of vision**—slender forehead support is unobtrusive and feels light on the face.
- **Designed for your lifestyle**—unique partition design in elbow effectively washes out carbon-dioxide and reduces noise ... so you and your partner can enjoy a good night's sleep.

## Intuitively easy

- **First-time fit**—it's easy to get acquainted with your mask when it fits you well the very first time you put it on. Your Mirage Micro is designed to fit ... first time, every time.
- **Ready to wear**—handy headgear clips are easy to attach and detach, and save your optimal headgear settings ... so you can set and forget!
- **On-the-face adjustment**—there's no need to take off the mask to use the MicroFit dial. You can simply adjust the dial with one hand while wearing the mask, and even while lying down.
- **Easy to care for your mask**—few and simplified parts make your mask simple to assemble, disassemble and clean.

## Mirage Micro

NASAL MASK

PROVEN SEAL, PERFECT FIT

