



Mirage™ SoftGel

NASAL MASK

Instant Softness, Lasting Comfort

The Mirage SoftGel Nasal Mask offers unprecedented softness and comfort. The DoubleGel cushion, with its clear gel layer provides comfort, while the blue gel layer maintains support and stability.

Backed by ResMed's renowned performance and quality, the frame accommodates two interchangeable cushions—Mirage SoftGel and Mirage Activa™ LT.

FOR THE PATIENT

Soft comfort meets stability

- DoubleGel cushion combines softness with stability.

A good night's sleep

- An array of small, built-in vents diffuses air quietly providing a peaceful environment for patients and their bed partners.

Fitting made easier

- The Mirage dual-wall with the DoubleGel cushion makes fitting simple.
- Cushion conforms to the unique facial contours of each individual, creating a comfortable fit.

Interchangeable cushions offer choice

- Choice of cushions makes it easier to remain compliant without completely replacing the mask.

FOR THE PROFESSIONAL

Customer satisfaction, added compliance

- Increased therapy acceptance through first-time fit—lasting comfort.
- Comprehensive user support package.

Improved efficiencies

- Interchangeable cushions and shared components improve inventory efficiencies, reducing costs, fitting and training times.
- A multi-patient, multi-use mask, Mirage SoftGel is suitable for repeated use in the lab.
- Extremely durable, Mirage SoftGel ensures consistent performance with fewer callbacks.

- **Soft, inviting, makes a great 'first impression'**
- **Excellent for new users**
- **Comfort throughout therapy**
- **Renowned quality and reliability**

PATIENT-PREFERRED

- 8 out of 10 patients rate Mirage SoftGel better in comfort and durability than other gel masks*.
- 3 out of 4 patients rate their sleep better with Mirage SoftGel*.

* ResMed User Evaluation Trials—Mirage SoftGel vs Competitor Gel Masks



MicroFit dial

offers 24 positions for better fit, seal and comfort

ResMed's unique frame

accommodates both Mirage SoftGel and Mirage Activa LT cushions for business efficiency, cost savings and flexibility

360° elbow rotation

positions the tubing for ease of use


Quiet, diffused vent design

disperses air softly and quietly, away from bed partner



Elbow partition

aids CO₂ washout and reduces noise output



Mirage SoftGel's unique DoubleGel cushion

Clear gel layer provides softness and comfort; blue gel layer provides support and stability

Set-and-forget headgear clips

retain optimal headgear settings; easy to attach and detach

Swivel

offers a convenient way to connect and disconnect air tubing

PRODUCT CODES

Mirage SoftGel Nasal Mask

Asia Pacific

Small	61642
Medium	61643
Large	61644
Large Wide	61645