



ResMed

**Actual size: 140mm (Width)
x 175mm (Height)**

Ruler is for scaling purposes.

**Use this fitting template as
a guide to size your patients.**

How to use the fitting template:

- Cut out the template along the dotted line
- Follow directions on template to determine the recommended cushion size.
- Please ensure the template is printed to scale.


Quattro™ Air

AirFit™ F10

Full face mask
Full face mask **for Her**

Rx Only 
0123


AirFit and Quattro are trademarks of ResMed Ltd. Quattro is registered in the U.S. Patent and Trademark Office.
© 2014 ResMed Ltd. 1017454/1 2014-04




Quattro™ Air

AirFit™ F10


Full face mask
Full face mask **for Her**




SIZE




AirFit and Quattro are trademarks of ResMed Ltd.
Quattro is registered in the U.S. Patent and Trademark Office.
© 2014 ResMed Ltd. All rights reserved. Specifications may change without notice.
Product availability may vary across regions.

Rx Only 
0123



ResMed



Sample size check