

RESMED



S7™ Lightweight

Patient Information

The value-added choice for any season


For looks and user-friendliness, the S7 Lightweight scores highly. With or without the humidifier, the unit is stylish and compact so it sits attractively on your bedside table.

Having a properly fitted mask that is comfortable and doesn't leak is vital for successful CPAP therapy. The S7 Lightweight is fully compatible with the complete range of ResMed masks. With a reputation for comfort and seal, ResMed masks are available in a range of models and sizes to suit all face shapes.

This stylish and compact unit sits attractively on your bedside table.



ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd North Ryde, NSW, Australia +61 2 9886 5000 or 1 800 658 189 (toll free). Offices in Austria, Finland, France, Germany, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, United Kingdom (see web site for details).

S7 Lightweight is protected by patents: US 4944310, US 5199424, US 5522382. Designs registered: AU 147283, AU 147335, AU 147336, CH 128.709, CH 128.710, CH 128.711, FR 02 1407, GB 3001791, GB 3001819, GB 3001820, GB 3001821, US D467335. Other designs pending. HumidAire 2i and  are trademarks of ResMed Ltd. © 2003 ResMed Ltd. 101978/1 03 03



a CPAP for all seasons



S7™ *Lightweight*

Introducing ResMed's new S7 Lightweight—the standard CPAP that gives you outstanding therapy and value. Delivering quiet, fixed, continuous positive airway pressure, the S7 Lightweight gives you the option to add integrated humidification—if and when you need it.

The flexible choice for humidification

The cool, dry air from CPAP therapy may make your nose, mouth, or throat dry, particularly in winter. It may also cause a runny nose, sneezing, or nasal congestion. Heated humidification can help relieve these symptoms by adding moisture to the air before you breathe it. As well as improving comfort, humidification can help to prevent mouth breathing if you have any of the above symptoms. Mouth breathing may decrease the effectiveness of your therapy, so it is important to manage this effectively.

The S7 Lightweight has an optional HumidAire 2i humidification module.

Small and portable, the HumidAire 2i simply clips onto the S7 Lightweight to give you truly integrated humidification—if and when you need it. It is ideal if you find you don't need humidification during the warmer months, or if you prefer not to take a humidifier when you travel. Your needs for

humidification can also vary from night to night.

The HumidAire 2i has a unique dial for adjusting the amount of humidification added to your CPAP pressure that is easy to use, even in the dark.

The reassuring choice for the first time user

Although humidification adds comfort and aids therapy, you may want to try CPAP therapy without humidification when you first start using it. The S7 Lightweight is a risk-free choice. You can simply purchase or rent the HumidAire 2i after experiencing therapy if you, or your physician, feel it is needed.

The convenient choice season after season

The S7 Lightweight is one of the lightest CPAP devices on the market, and has a convenient carry handle that allows the unit to be carried with just one hand.

Traveling with the S7 Lightweight is also made easy thanks to a universal mains supply that automatically switches to local power. Altitude compensation is included for more accurate pressure delivery when traveling, and the device comes with a tough travel bag for protection.



Simply clip on the HumidAire 2i—if and when you need it

Easy to use modules make the S7 Lightweight truly flexible