



# Data management/Device compatibility list

English

This supplement specifies the devices compatible with ResMed data management products. For full details regarding the correct use of these devices, see the relevant User Guides.

	ResScan™ 5.9	EasyCare Tx™ / Tx Link™ 8.6	AirView™	RCM	RCMH**
<b>AirMini™</b>			✓		
<b>Lumis Series</b>					
Lumis™ 100 VPAP S	✓	✓	✓		
Lumis™ 100 VPAP ST	✓	✓			
Lumis™ 150 VPAP ST	✓	✓	✓		
Lumis™ 100 VPAP ST-A	✓	✓	✓		
Lumis™ 150 VPAP ST-A	✓	✓	✓		
<b>AirSense 10 / AirCurve 10 Series*</b>					
AirSense™ 10 AutoSet	✓	✓	✓		
AirSense™ 10 AutoSet for Her	✓	✓	✓		
AirSense™ 10 Elite	✓	✓	✓		
AirSense™ 10 CPAP	✓		✓		
AirSense™ 10 Auto	✓	✓	✓		
AirSense™ 10 Respond A	✓	✓	✓		
AirSense™ 10 Respond W	✓	✓	✓		
Sleepmate 10 Auto	✓	✓	✓		
Sleepmate 10 Elite	✓	✓	✓		
AirCurve™ 10 CS PaceWave™			✓		
<b>S9 Series</b>					
Auto 25 (S9)	✓	✓			
AutoSet CS-A (S9 PaceWave™)	✓	✓			
AutoSet CS™ (S9 PaceWave™)	✓	✓			
AutoSet CS™ (S9)	✓	✓			
S9 AutoSet™	✓	✓	✓		
S9 Elite™	✓	✓	✓		
S9 Escape™	✓	✓	✓		
S9 Escape™ Auto	✓	✓			
S9 VPAP™ Tx		✓			
V-Auto (S9)	✓	✓			
VPAP™ Adapt (S9 PaceWave™)	✓	✓			
VPAP™ Adapt (S9)	✓	✓			
VPAP™ ST (S9 iVAPS)	✓	✓			
VPAP™ ST-A (S9 iVAPS)	✓	✓			
VPAP™ ST-A (S9)	✓	✓			
VPAP™ S (S9)	✓	✓			
VPAP™ ST (S9)	✓	✓			
Lumis Tx	✓				
RESMED S9 Auto	✓				
S9 Respond	✓				
Sleepmate S9 Auto	✓				
Sleepmate S9 Elite	✓				



	ResScan™ 5.9	EasyCare Tx™ / Tx Link™ 8.6	AirView™	RCM	RCMH**
<b>S8 Series</b>					
S8 AutoScore™	✓	✓			
S8 AutoScore™ II	✓	✓			
S8 AutoSet Spirit™	✓	✓			
S8 AutoSet Spirit™ II	✓	✓			
S8 AutoSet Vantage™	✓	✓			
S8 AutoSet™ II	✓	✓			
S8 Elite™	✓				
S8 Elite™ II	✓	✓			
S8™ Auto 25	✓	✓			
VPAP™ IV Series	✓	✓			
<b>S7 Series</b>					
AutoSet CS™ 2	✓	✓			
S7 AutoSet Spirit™	✓				
S7 Elite™	✓				
VPAP™ Adapt SV	✓				
VPAP™ II Series	✓				
VPAP™ III Series	✓	✓			
VPAP™ III ST-A with QuickNav	✓	✓			
<b>Ventilators</b>					
Astral™ 100	✓			✓	✓
Astral™ 150	✓			✓	✓
Stellar™ 100	✓	✓		✓	
Stellar™ 150	✓	✓		✓	
Elisée™ 150	✓				
<b>Home sleep testing device</b>					
ApneaLink™ Air			✓		

\* Compatible with AirView via wireless connectivity only.

\*\* For a list of available data for use by the EMR system, please contact your ResMed representative.

RCM – ResMed Connectivity Module

RCMH – ResMed Connectivity Module Hospital

Notes:

- Not all products available in all regions.
- AirView compatibility may not be available in all regions.

 **ResMed Ltd** 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirMini, AirSense, AirView, Astral, ApneaLink, Lumis, S8, S9, Stellar, VPAP, AutoSet, S9 Escape, S9 Elite, S9 Escape Auto, S8 Elite, S8 AutoSet Vantage, PaceWave, ResScan, VPAP Tx, S8 AutoSet, S9 AutoSet, S9 Auto 25, VS III, Elisée, VPAP Adapt SV, AutoSet CS, S8 AutoScore, S8 AutoSet Spirit and VPAP are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip. © 2019 ResMed Ltd. 1018136/13 2019-01