



# New Maxi-Therm® & Maxi-Therm® Lite Blanket Connectors

CSZ's Maxi-Therm® & Maxi-Therm® Lite blankets now have self-sealing connectors. No longer is there a need to "clamp" off the tubing of the blanket before disconnecting it from the device.



The new self-sealing connections help prevent water from coming out of the blanket or dripping from the clamps. CSZ believes this enhancement will make your job easier and your patients more comfortable. Nothing else has changed on our Maxi-Therm® or Maxi-Therm® Lite blanket, both still offer the same great therapy for your patient as before.

It is important to ensure that all of the right parts of your hyper-hypothermia system are in place. It is the integrity of the system (water blankets, hoses and device) that are most important when it comes to heat exchange to the patient. CSZ designs all of its products to work as a system and provides the proper therapy that your patient needs.

Cat. Number	Product Name	Description
870	Maxi-Therm® Lite	25" x 4" (31.8cm x 45.7cm)
871	Maxi-Therm® Lite	22" x 17.25" (55.9cm x 43.82cm)
872	Maxi-Therm® Lite	25" x 19" (63.5cm x 31.8cm)
873	Maxi-Therm® Lite	12.5" x 18" (31.8cm x 45.7cm)
874	Maxi-Therm® Lite	25" x 33" (63.5cm x 83.8cm)
876	Maxi-Therm® Lite	25" x 64" (63.5cm x 162.6cm)
273	Maxi-Therm®	12" x 18" (30.5cm x 45.7cm)
274	Maxi-Therm®	22" x 30" (55.9cm x 83.8cm)
276	Maxi-Therm®	24" x 60" (61cm x 152.4cm)