



Gyrolab Loading List

Run Setup name ADA Exp 1 - Initial Screen BA200
Created from version 1

Gyrolab Bioaffy 200



MP-1



Gyrolab Loading List

Microplate	Type
MP-1	MP96 PCR

Pos	Reagent/Sample	Conc.	Spike conc.	Dilution	Vol. [μ l]	Well
1	PC MM3_without drug_Blank	0,0 ng/mL		1,0	4	A1
2	PC MM3_without drug_1	400,0 ng/mL		1,0	4	A2
3	PC MM3_without drug_2	200,0 ng/mL		1,0	4	A3
4	PC MM3_without drug_3	100,0 ng/mL		1,0	4	A4
5	PC MM3_without drug_4	50,0 ng/mL		1,0	4	A5
6	PC MM3_without drug_5	25,0 ng/mL		1,0	4	A6
7	PC MM8_with drug_Blank	0,0 ng/mL		1,0	6	B1
8	PC MM8_with drug_1	400,0 ng/mL		1,0	4	B2
9	PC MM8_with drug_2	200,0 ng/mL		1,0	4	B3
10	PC MM8_with drug_3	100,0 ng/mL		1,0	4	B4
11	PC MM8_with drug_4	50,0 ng/mL		1,0	4	B5
12	PC MM8_with drug_5	25,0 ng/mL		1,0	4	B6
13	PC MM3_with drug_Blank	0,0 ng/mL		1,0	6	C1
14	PC MM3_with drug_1	400,0 ng/mL		1,0	4	C2
15	PC MM3_with drug_2	200,0 ng/mL		1,0	4	C3
16	PC MM3_with drug_3	100,0 ng/mL		1,0	4	C4
17	PC MM3_with drug_4	50,0 ng/mL		1,0	4	C5
18	PC MM3_with drug_5	25,0 ng/mL		1,0	4	C6
19	MM3 No Drug			1,0	4	D1
20	MM3 Conf. Drug level 1			1,0	4	D2
21	MM3 Conf. Drug level 2			1,0	4	D3
22	MM8 No Drug			1,0	4	D4
23	MM8 Conf. Drug level 1			1,0	4	D5
24	MM8 Conf. Drug level 2			1,0	4	D6
25	Wash Solution				100	E1
26	Wash Solution				100	E2
27	Wash Solution				100	E3

Gyrolab Loading List

MP-1

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	2	3	4	5	6						
B	7	8	9	10	11	12						
C	13	14	15	16	17	18						
D	19	20	21	22	23	24						
E	25	26	27									
F												
G												
H												