


Gyrolab Loading List

Run Setup name ADA Exp 3 - Cut point Mix96
Created from version 1

Gyrolab Mixing CD 96



MP-1

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	2	3	4	5	6	7	8	9	10	11	12
B	13	14	15	16	17	18	19	20	21	22	23	24
C	25	26	27	28	29	30	31	32	33	34	35	36
D	37	38	39	40	41	42	43	44	45	46	47	48
E	49	50										
F	51	52	53	54								
G												
H												

	Blank		QC		Unknown		Reagent
	Standard		NC/PC				Wash solution
	CC		Affinity				

Gyrolab Loading List

Microplate	Type
MP-1	MP96 PCR

Pos	Reagent/Sample	Conc.	Spike conc.	Dilution	Vol. [μ l]	Well
1	PC 1 __ Spiked			1,0	5	A1
2	PC 1 __ Unspiked			1,0	5	A2
3	PC 2 __ Spiked			1,0	5	A3
4	PC 2 __ Unspiked			1,0	5	A4
5	PC 3 __ Spiked			1,0	5	A5
6	PC 3 __ Unspiked			1,0	5	A6
7	NC 1 __ Spiked			1,0	5	A7
8	NC 1 __ Unspiked			1,0	5	A8
9	NC 2 __ Spiked			1,0	5	A9
10	NC 2 __ Unspiked			1,0	5	A10
11	NC 3 __ Spiked			1,0	5	A11
12	NC 3 __ Unspiked			1,0	5	A12
13	Sample1_Sample1__ Spiked				5	B1
14	Sample1_Sample1__ Unspiked				5	B2
15	Sample2_Sample2__ Spiked				5	B3
16	Sample2_Sample2__ Unspiked				5	B4
17	Sample3_Sample3__ Spiked				5	B5
18	Sample3_Sample3__ Unspiked				5	B6
19	Sample4_Sample4__ Spiked				5	B7
20	Sample4_Sample4__ Unspiked				5	B8
21	Sample5_Sample5__ Spiked				5	B9
22	Sample5_Sample5__ Unspiked				5	B10
23	Sample6_Sample6__ Spiked				5	B11
24	Sample6_Sample6__ Unspiked				5	B12
25	Sample7_Sample7__ Spiked				5	C1
26	Sample7_Sample7__ Unspiked				5	C2
27	Sample8_Sample8__ Spiked				5	C3
28	Sample8_Sample8__ Unspiked				5	C4
29	Sample9_Sample9__ Spiked				5	C5
30	Sample9_Sample9__ Unspiked				5	C6
31	Sample10_Sample10__ Spiked				5	C7
32	Sample10_Sample10__ Unspiked				5	C8
33	Sample11_Sample11__ Spiked				5	C9
34	Sample11_Sample11__ Unspiked				5	C10
35	Sample12_Sample12__ Spiked				5	C11
36	Sample12_Sample12__ Unspiked				5	C12
37	Sample13_Sample13__ Spiked				5	D1
38	Sample13_Sample13__ Unspiked				5	D2
39	Sample14_Sample14__ Spiked				5	D3
40	Sample14_Sample14__ Unspiked				5	D4
41	Sample15_Sample15__ Spiked				5	D5
42	Sample15_Sample15__ Unspiked				5	D6
43	Sample16_Sample16__ Spiked				5	D7
44	Sample16_Sample16__ Unspiked				5	D8
45	Sample17_Sample17__ Spiked				5	D9

Gyrolab Loading List

Pos	Reagent/Sample	Conc.	Spike conc.	Dilution	Vol. [μ l]	Well
46	Sample17_Sample17__Unspiked				5	D10
47	Sample18_Sample18__Spiked				5	D11
48	Sample18_Sample18__Unspiked				5	D12
49	Acidic Buffer_0.5 M Glycine-HCl pH 2.6				85	E1
50	Master Mix_Neutralization Buffer				85	E2
51	Wash Solution				100	F1
52	Wash Solution				100	F2
53	Wash Solution				100	F3
54	Wash Solution				100	F4

MP-1

