



Gyrolab Loading List

Run Setup name ADA Exp 4 - Performance BA200
Created from version 1

Gyrolab Bioaffy 200



MP-1



Gyrolab Loading List

Microplate	Type
MP-1	MP96 PCR

Pos	Reagent/Sample	Conc.	Spike conc.	Dilution	Vol. [μ l]	Well
1	Drug dil - without PC- MM X_Diluted pool serum_Blank	0,0 ug/mL		1,0	5	A1
2	Drug dil - without PC- MM X_Diluted pool serum_1	1280,0 ug/mL		1,0	4	A2
3	Drug dil - without PC- MM X_Diluted pool serum_2	640,0 ug/mL		1,0	4	A3
4	Drug dil - without PC- MM X_Diluted pool serum_3	320,0 ug/mL		1,0	4	A4
5	Drug dil - without PC- MM X_Diluted pool serum_4	160,0 ug/mL		1,0	4	A5
6	Drug dil - without PC- MM X_Diluted pool serum_5	80,0 ug/mL		1,0	4	A6
7	Drug dil - without PC- MM X_Diluted pool serum_6	40,0 ug/mL		1,0	4	A7
8	Drug dil - without PC- MM X_Diluted pool serum_7	20,0 ug/mL		1,0	4	A8
9	Drug dil - without PC- MM X_Diluted pool serum_8	10,0 ug/mL		1,0	4	A9
10	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_Blank	0,0 ug/mL		1,0	5	B1
11	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_1	1280,0 ug/mL		1,0	4	B2
12	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_2	640,0 ug/mL		1,0	4	B3
13	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_3	320,0 ug/mL		1,0	4	B4
14	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_4	160,0 ug/mL		1,0	4	B5
15	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_5	80,0 ug/mL		1,0	4	B6
16	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_6	40,0 ug/mL		1,0	4	B7
17	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_7	20,0 ug/mL		1,0	4	B8
18	Drug dil - 100 ng/mL PC - MM X_dil pool serum + 100 ng/mL PC_8	10,0 ug/mL		1,0	4	B9

Gyrolab Loading List

Pos	Reagent/Sample	Conc.	Spike conc.	Dilution	Vol. [μ l]	Well
19	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_Blank	0,0 ug/mL		1,0	5	C1
20	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_1	1280,0 ug/mL		1,0	4	C2
21	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_2	640,0 ug/mL		1,0	4	C3
22	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_3	320,0 ug/mL		1,0	4	C4
23	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_4	160,0 ug/mL		1,0	4	C5
24	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_5	80,0 ug/mL		1,0	4	C6
25	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_6	40,0 ug/mL		1,0	4	C7
26	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_7	20,0 ug/mL		1,0	4	C8
27	Drug dil - 800 ng/mL PC - MM X_dil pool serum + 800 ng/mL PC_8	10,0 ug/mL		1,0	4	C9
28	PC without drug - MM X_Diluted pool serum_Blank	0,0 ng/mL		1,0	6	D1
29	PC without drug - MM X_Diluted pool serum_1	800,0 ng/mL		1,0	4	D2
30	PC without drug - MM X_Diluted pool serum_2	400,0 ng/mL		1,0	4	D3
31	PC without drug - MM X_Diluted pool serum_3	200,0 ng/mL		1,0	4	D4
32	PC without drug - MM X_Diluted pool serum_4	100,0 ng/mL		1,0	4	D5
33	PC without drug - MM X_Diluted pool serum_5	50,0 ng/mL		1,0	4	D6
34	PC without drug - MM X_Diluted pool serum_6	25,0 ng/mL		1,0	4	D7
35	PC without drug - MM X_Diluted pool serum_7	12,5 ng/mL		1,0	4	D8
36	Wash Solution				100	E1
37	Wash Solution				100	E2
38	Wash Solution				100	E3
39	Wash Solution				100	E4

Gyrolab Loading List

MP-1

	1	2	3	4	5	6	7	8	9	10	11	12
A	1	2	3	4	5	6	7	8	9			
B	10	11	12	13	14	15	16	17	18			
C	19	20	21	22	23	24	25	26	27			
D	28	29	30	31	32	33	34	35				
E	36	37	38	39								
F												
G												
H												