KMO[™]V-V Knee Orthosis

Features

- No need to contour uprights with Single posterior joint design
- Patient friendly terry cloth liner system
- Pre and post operative management of bony and soft tissue injuries
- Adjustable range of stabilization of the knee in the sagittal and/or the coronal plane
- Static progressive positioning from 0° to 110° (additional ROM may be attained depending on size of patient's leg)
- Adjustable length from 12.50" to 16.50"(Pediatric) and from 17" to 26"(Adult)
- Available in Adult and Pediatric models, along with additional KMO liner variations

Other KMO[™]V-V Variations



KMO[™]V-V (Knee Management Orthosis with Varus Valgus adjustment), with it's single posterior joint, is a cost-effective, patient friendly, and a simple approach for static progressive positioning management of the knee.



KMO™ (w/ PLASTIC THIGH & CALF SHELL) KNEE ORTHOSIS

Uses

- Chronic systemic conditions
- » Immediate prophylactic protection at end range until treatment plans are finalized
- Contracture management
- ▶ CVA, TBI, SCI, CP, and MS
- Emergent lower limb immobilization
- Musculoskeletal injuries
- BKA(Below the Knee Amputation)
- Post Botox or serial cast management
- TKA(Total Knee Arthroplasty)
- Post-operative immobilization

PDAC APPROVED L-1831



