

EV™ (Equino-Varus) Ankle Foot Orthosis

Positioning and Walking Adjustable Hinged Ankle Foot
Brace with Inversion and Eversion Adjustable Joint

Features-

- › Offers positive heel suspension
- › Accommodates and helps control varus and valgus conditions of the ankle/foot complex
- › The heel bar features calibrated settings for a more accurate documentation of the patients progress
- › Has an adjustable anti-rotation bar which helps control unwanted extremity rotation
- › The EV™ Orthosis bar works very well in standard AFO's, but also is extremely useful when used as the distal segment (PKA™ Orthosis) in a custom-made KAFO or a pre-fabricated KO or KAFO system.



*** Heel strap positioned inside the heel connecting bar for display purposes only*



Uses-

- › With its securely fastened walking base, the EV™ Articulating Ankle/Foot Orthosis can be used on the ambulatory as well as the recumbent patient.



****Available in Adult #651SKT, Pediatric #551SKT, and Infant #451RKG**