

It's more than just positioning and offloading!

RIBBY™

The RIBBY™ Orthosis is a fantastic and creative way to display your facility's ability to provide excellent orthotic care that addresses the rehab patient's needs 24/7, in a most cost-effective and yet profitable manner.

Features

- Clinically proven for effective ambulation*
 - Dorsiflexion adjustability (+ or - 20°) Plantarflexion adjustability (+ or - 40°)
 - Quick disconnect free dorsiflexion/plantarflexion stop-joint
 - Moldable foot plate
 - Open toe shoe provides at no extra cost
 - Compatible with all liner variations that we offer if applicable Proven superior integrity of the superstructure**
- Available in Adult Standard # 656NPB



Uses

- Rehab patients in orthotic care.
- Early intervention for the evolving rehab patient
 - CVA, TBI, SCI and Diabetic neuropathy
- Facilitation of functional knee and ankle joint position
 - Allows for optimal dorsi/plantar flexion ankle positioning to provide adequate clearance during swing phase and knee stability in stance phase.
- Post musculoskeletal injury of the foot and ankle
 - Post TAL immobilization
 - Post Botox or serial casting
 - Stable fractures of the foot and ankle
 - Tendon or ligamentous tear/rupture

**DeToro W, (2001) Plantarflexion Resistance of Selected Ankle-Foot Orthoses: A Study of Commonly Prescribed Prefabricated and Custom-Molded Alternatives: J Prosthet Orthot:13;39-44

PDAC Approved L-1971

