



caleva

PRODUCT INFORMATION

LABORATORY

CALEVA

MIXER - GRANULATOR - BLENDER

Effectively mix, blend or granulate small batches of your powders,
wet masses, gels, waxes, slurries and pastes etc. on your bench top

PRODUCT INFORMATION

One unit, many uses



CALEVA CALEVA MIXER, GRANULATOR, BLENDER

- ✓ Imagine an effective mixer, blender and granulator for very small quantities on your laboratory bench-top.
- ✓ There is no other laboratory equipment that is easy to work with whilst showing an equivalent versatility and effectiveness with small batches.
- ✓ A range of different mixing blade / bowl configurations and inserts are available to allow the Caleva MGB to work with the widest possible range of products and "difficult to mix" materials.
- ✓ Mixing and blending is rapid in operation making your development faster and more likely to succeed.
- ✓ Save money by working efficiently with mixing, blending and granulation of your small batches with a wide variety of materials



Two sizes of bowl are available to suit your requirement

DESCRIPTION

A component of the Caleva Multi Lab that has found its own place as a separate and versatile small laboratory Mixer/Blender/Granulator, is the Caleva MGB. It is the best choice for a laboratory that needs to complete these processes quickly with small sample sizes. An ideal mixer for many development processes designed to improve the chances of success in development projects.


The small size and ease of cleaning reduces the time to complete trial iterations, offering faster results to enable you to reach your objectives more quickly.
Speed up your R&D, make the best use of your available budgets.

Different mixing bowl and blade sets as well as inserts and water jackets are available, and the best choice will depend on the properties of your product. Test your products with us to assist in your decision making. [Contact us to discuss what you want to achieve.](#)

Water jacketed bowls can be used to add or remove heat from the mixing system.

Different bowl sizes are available for a range of batch sizes. Bowl inserts can be used to further enhance your mixing options. Talk to us and let us help you make the right choices.

Test your product before you decide which system components are right for you so that the Caleva MGB will do what you want it to do.



A wide range of bowl and blade configurations are available depending on your requirements

WILL IT HELP YOUR LABORATORY TO ACHIEVE YOUR OBJECTIVES?

We have:


- Tested, and been able to mix effectively, a wide range of products from different industries from food and pharmaceutical to catalysts and petrochemicals as well as inorganic chemicals and crystalline salts.
- Blended and distributed a wax at room temperature evenly distributing this into a light natural organic powder.
- Blended and mixed a wide range of different pastes and gels.
- The list goes on and we are still finding new applications and industries for this flexible laboratory unit.

No matter what your industry or your product characteristics talk to us, discuss what you want to achieve and let's see if this laboratory mixer/blender/granulator (the Caleva MGB) will offer you the additional capability you need in your laboratory.

Want to test your products? [Contact us!](#)

THE CALEVA MGB IN USE

- Provides your laboratory with an opportunity to mix, granulate and blend a wide range of small samples.
- Development progress can be faster and new ideas can be tested quickly with minimum use of product, increasing chances of R&D success.
- Complete more iterations of your development trials and achieve your objectives faster.



Bowl inserts are available for highly efficient mixing of very small batch sizes

- A bench-top mixer, blender and granulator with lots of options all in a single, easy to use unit.
- ✓ **Make it yours!**
 - Have more space on your bench-top
 - Optimize your R&D
 - Increase your productivity
 - Complete more trials with fast and easy cleaning

TALK TO US

Please call us without obligation

+44 (0) 1258 471122

info@caleva.com



Certificate Number: 1503