


# 12 REASONS SCAN USERS WILL LOVE THE LANTOS SYSTEM



1

## SCAN SURGICAL PATIENTS

Surgical ears? No problem. Our membrane-based scanning system eliminates the need to pack surgical cavities.

2

## CAPTURE DEEP EAR CANALS

The Lantos membrane sits just millimeters from the eardrum, so you can capture deep ear canals for IICs easily and with no patient discomfort.

3

## SEE WHERE YOU'RE GOING

With full visualization into the ear canal, you (and your patients) can see where you are at every step in the ear scanning process.

4

## NO MORE OTOBLOCKS

Our membrane-based system doesn't use otoblocks, so no more worries about otoblock insertion or removal.

5

## NO RISK OF BLOW-BY

Without the use of silicone, there's no risk of impression material ending up where it shouldn't.

6

## REDUCE THE NEED FOR WAX REMOVAL

As the membrane inflates with a water-based solution, it gently pushes aside hair and small amounts of wax. This means no problem scanning ears that contain non-occluding wax.

7

## RETAIN THE SCAN ON FILE

Did a customer lose a hearing aid? Do they need sleep plugs or swim plugs? With a digital scan on file, orders and re-orders are effortless, with no need for additional patient visits.

8

## ONE SCAN, ENDLESS OPTIONS

Whether you're planning to prepare a single custom-fit hearing device or the full array of hearing enhancement and protection solutions, a single scan can be used for unlimited orders.

9

## DIGITAL IMPRESSION STORAGE

Ditch that cabinet of physical impressions that you keep "just in case"!

10

## NO MORE SHIPPING OF EMI

The digital ear scan gets uploaded electronically for manufacturing, so no more packaging, sorting, or shipping physical earmold impressions.

11

## SPEED UP PRODUCT TURNAROUND

By eliminating the need to ship physical earmold impressions, the Lantos System means your patients get their custom-fit products 2-4 days faster.

12

## MORE COMFORTABLE FOR PATIENTS

Not only is the Lantos 3D Ear Scanning more efficient than silicone-based earmold impressions, it is also much more comfortable for patients. In fact, user surveys indicate that more than 9 out of 10 patients prefer ear scanning to EMI.

## ★ BONUS: GET CEUs ★

Complete the Lantos training and receive up to .6 CEUs.

LANTOS  
TECHNOLOGIES

