

EndoSheath® Technology

The sterile, disposable solution for flexible endoscopy

EndoSheath® Technology is not your typical barrier. Designed to offer a strong, durable, proven effective microbial barrier, this is not just a simple cover. By providing a barrier between the endoscope and patient, EndoSheath® Technology allows for less complicated cleaning and disinfection routines, and improves equipment turnaround time. All the difficult-to-clean endoscope components with EndoSheath® Technology are disposable, including: the barrier, the channels, the ports, and the seals. Clinically proven, EndoSheath® Technology allows for efficient, effective endoscope reprocessing in any setting.

With EndoSheath® Technology there is strength in numbers...


- › 25+ years on the market
- › 30+ FDA clearances
- › 37+ systems with CE Mark
- › 9+ systems with HealthCanada license
- › 5+ million procedures performed worldwide
- › 0 cross-contamination complaints



Reusable Endoscope



with sterile, disposable
EndoSheath® Technology



EndoSheath®
ENDOSCOPY SYSTEMS

A brand with strength. A name you can trust.

Now included
in the
2015 AAMI/ANSI
National Standards
for Reprocessing

Pioneering Infection Control with EndoSheath® Endoscopy

EndoSheath® Advantages:

- EndoSheath® Technology is a proven effective barrier to organisms as small as 27 nanometers per FDA requirements and testing
- Over 5 Million EndoSheath® products sold over 25 years without a single reported complaint of patient-to-patient cross contamination
- EndoSheath® Technology is designed to stretch and maintain its integrity, and not to tear or break during procedures
- All EndoSheath® disposables are 100% leak tested during manufacturing with a 0% acceptable failure rate for production
- All barriers, channels, ports, seals, and tubing are sterile and disposable
- Concerns over biofilm, bioburden, and issues of improperly rinsed devices are eliminated for patient contact areas, as they are always new, sterile, and disposable
- Easy to use and prepare limiting the complexity of endoscopy preparation routines
- Compatible with standard accessories such as biopsy forceps, needles, and other endoscopic devices
- Over 30 research studies performed on EndoSheath® Technology since 1992
- Disposables feature a 3-year shelf life
- Included in the ANSI/AAMI ST91:2015 National Standards for Flexible and semi-rigid endoscope processing in health care facilities
- EndoSheath® Technology is FDA Cleared, CE Marked, and CDC Recognized

Currently available for:
Bronchoscopy, Cystoscopy,
Esophagoscopy, and Laryngoscopy

Cogentix Medical, Inc.

40 Ramland Road South / Orangeburg, NY 10962

Telephones: 800-874-9975 / 845-365-0600 / Fax: 845-365-0620

E-Mail: info@cogentixmedical.com



www.cogentixmedical.com

Innovative Solutions:

Sterile, durable, disposable
Microbial barrier




A reusable endoscope

A sterile, disposable channel
for every procedure



Unique "D-Shape"
endoscopes



Patented system allows
for complete isolation of
endoscope from patient

About Cogentix Medical, Inc.

Cogentix Medical, Inc. has brought pioneering concepts to the world of flexible endoscopy for over 20 years. With the patented EndoSheath® Technology and innovative endoscope designs, Cogentix Medical, Inc. has transformed the economics of flexible endoscopy while enhancing patient safety for Otolaryngology, Gastroenterology, Urology, and Pulmonology practitioners. Cogentix Medical has made a commitment to help physicians transform their practices with the EndoSheath® Technology.

Request a **Clinical Research Summary** or an Official EndoSheath® Technology **Infection Control package:**

Call **800-874-9975** or visit **www.cogentixmedical.com**