



Single-Type AV Fistula Needle

JMS introduces Single-Type AV Fistula needles in Sharp and Buttonhole styles.




JMS Single-Type AV Fistula needles are offered in 20cm and 40cm tube length. The longer sets are designed for self-cannulators to easily access upper arm fistulas. The Harmony[®] buttonhole sets include our unique site preparation tool.


Standard Needle



Low recirculation volume in both styles.



Back-eye needle provides for efficient flow rates.




Sharp set includes familiar Red and Blue clamps.


Harmony®



Includes site preparation tool intergrated into needle cover.



Low recirculation volume and back-eye needle.



Yellow clamps indicate buttonhole needle.

***JMS Single-Type requires HD equipment capable of Single Lumen, Single Needle therapy**

Item Code	Description	Case Pack
820-3500-31	15Gx1" AVF Single-Type, Red/Blue Clamps, Sharp	500
820-3600-31	16Gx1" AVF Single-Type, Red/Blue Clamps, Sharp	500
820-3507-31	15Gx1" AVF Single-Type, Red/Blue Clamps, Sharp, 40cm	250
820-3605-31	16Gx1" AVF Single-Type, Red/Blue Clamps, Sharp, 40cm	250
827-1507-31	15Gx1" AVF Single-Type, Yellow Clamps, Button-Hole, 40cm	250
827-1605-31	16Gx1" AVF Single-Type, Yellow Clamps, Button-Hole, 40cm	250



For more information please contact us or your local distributor

Copyright © 2016 JMS North America Corporation
All rights reserved

BRO/STAVF-E/001

**Imported By:
JMS North America Corporation**

MAIN OFFICE

22320 Foothill Blvd., Suite 350
Hayward, CA 94541, U.S.A.
Tel: +1-510-888-9090
Fax: +1-510-888-9099
NA: customerservice@jmsna.net

SALES & MARKETING

480 Sawgrass Corporate Pkwy, Suite 120
Sunrise, FL 33325, U.S.A.
Tel: +1-954-689-9280
Fax: +1-954-689-9289
LatAm: ventas@jmsna.net