

SysLoc[®]MINI

JMS introduces the **NEW** SysLoc[®]MINI

The latest design in safety AVF Needles. The SysLoc[®]MINI combines innovative technology and the audible "click" to ensure safety for patients and health-care workers.



SysLoc[®]MINI is over 20% shorter than the previous model providing easier use.



Shorter transparent hub improves control while providing full view of flash-back for proper positioning.



Contoured rotating hub provides for optimal flow and patient comfort.




The audible click ensures that the needle is securely locked into place – reducing the chance for needle-stick injury.



www.jmsna.net

JMS

System Safety Activation



For additional operating instructions please refer to the instructions for use.

Ordering Information

Single Set

Item Code	Description	Each / Case
864-1400-33	AVF Set 14G X 1" BE w/Clamp 12" Tubing	300
864-1500-33	AVF Set 15G X 1" BE w/Clamp 12" Tubing	300
864-1600-33	AVF Set 16G X 1" BE w/Clamp 12" Tubing	300
864-1700-33	AVF Set 17G X 1" BE w/Clamp 12" Tubing	300
864-1800-33	AVF Set 18G X 1" BE w/Clamp 12" Tubing	300

Twin Set

Item Code	Description	Twin Set / Case
864-2400-33	AVF Twin Set 14G X 1" BE w/Clamp 12" Tubing	250
864-2500-33	AVF Twin Set 15G X 1" BE w/Clamp 12" Tubing	250
864-2600-33	AVF Twin Set 16G X 1" BE w/Clamp 12" Tubing	250
864-2700-33	AVF Twin Set 17G X 1" BE w/Clamp 12" Tubing	250
864-2800-33	AVF Twin Set 15G X 1" BE w/Clamp 12" Tubing	250

*Other configurations are available upon request

For more information please contact us or your local distributor

Copyright ©2007 JMS North America Corporation.
All rights reserved.
BRO/SLM/001_R3

Imported By: JMS North America Corporation

Main Office
22320 Foothill Blvd., Suite 350
Hayward, CA 94541
Tel: 510-888-9090
Fax: 510-888-9099
NA: customerservice@jmsna.net

Sales & Marketing
480 Sawgrass Corporate Pkwy., Suite 120
Sunrise, FL 33325
Tel: 954-689-9280
Fax: 954-689-9289
LatAm: ventas@jmsna.net