

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: anticipating trends to help improve the quality of life.



4200 SERIES PNEUMATIC STOOLS

4200

SERIES PNEUMATIC STOOLS

Finished for Durability

- Select from 14 easy-to-maintain vinyl colours. Other vinyl colours available by request.
- All exposed metal, with exception of the aluminum base, is chrome-plated steel.
- Base and legs finished in tough, black epoxy enamel.
- Hooded dual-wheel casters.
- 5 legs for stability.

Optional Foot Ring

Order a foot ring with your stool, or at a later date.

Tech Specs

Base diameter: 22.5"



	4240	4246
Seat adjustment	Pneumatic Assist	Pneumatic Assist
Seat diameter	38 cm	38 cm
Seat construction	5 cm polyfoam	5 cm polyfoam
Seat swivels	360°	360°
Five legs for extra stability	Yes	Yes
Height range	49.5 - 68.6 cm	49.5 - 68.6 cm
Shipping weight	10.4 kg	11.8 kg

RA21062019/CP&BR

 HAAG-STREIT
INTERNATIONAL

 **DEVICE
TECHNOLOGIES**

📞 1300 DEVICE (338 423)

✉ customers@device.com.au

🌐 www.device.com.au