

Frey

SLI-200 Imaging System



SLI-200 Imaging System

SLI-200 Digital Imaging System - Visualise. Improve workflow. Share.

Frey SLI-200 provides great quality with a fully integrated design at the manufacturing stage and does not require an external PC. Implement efficient and high-level work flow into clinical practice with a less complicated, more integrated device.

Stream to multiple screens

Stream images and videos across multiple devices. Frey SLI-200 allows you to collaborate with other healthcare professionals, share knowledge, discuss diagnosis and improve clinical outcomes in real time, easily. Designed for general practices, educational and training institutions. Video streaming and photo sharing enable clinical collaboration and education.

Stunning AMOLED tech

The AMOLED screen provides superb image quality, achieved by its outstanding contrast and brightness as well as its impressive depth of colour. In-addition the wide viewing angle makes it ideal for viewing detailed images of the eye, whilst being positioned at the slit lamp.

Work wirelessly

Use fast, stable and reliable Wi-Fi AC or Bluetooth 4.2 to connect all your devices. Simplify workspace and improve clinical workflows.

HIGHLIGHTS	
Device type	SLI-200 Imaging System
Device features	AMOLED 5' touch screen
	Video streaming to multiple devices
	CMOS Image Sensors with a 3.45 µm Pixel
	128 GB internal high endurance flash storage
	HDR & Global Shutter Function for high-picture-quality imaging
	Photo & video capture
	Photos: 2048x1536
	HD video: 960x720
	Bluetooth 4.2/ Wi-Fi AC
	Colour Depth: up to 12 BPP
	Automatic L/R eye detection

Out of the box & stand-alone

Work independently and efficiently. Frey SLI-200 does not require an external PC or any other device to operate. All cables remain hidden in an ergonomic and minimalistic design. The eye lighting system, which can be controlled using the knob built into the lamp base, is included.

Large & reliable internal storage

Result and image management is an essential part of a doctor's workflow. SLI-200 simplifies and speeds up the work with 128GB of internal storage. Gather a large number of heavy files, select the most interesting photos and send them to your smartphone without cluttering it.

