

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



520

STYLISH ASSISTED RECLINE TILT CHAIR

520

ASSISTED RECLINE TILT EXAM CHAIR



Ergonomic Features

- Detached footswitch for hands-free raising/lowering.
- Large padded headrest keeps patients comfortable while reclined.
- Upholstered armrests rotate up to allow patients to move in and out from either chair side.
- Dual-rotation locks allow you to rotate chair 360° from a single position.
- Large hinged footrest supports patient's full weight, and folds away to save space.

Effortless Recline to 40°

- A manual recline bar, easily accessible from both sides of the chair, controls the assisted tilt of the chair up to 40°.

Power Hydraulic Base Boosts Efficiency

- The 520 chair's corded footswitch allows you to raise and lower the patient while maintaining contact.

Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam that lasts for years.
- Removable cushions for easy cleaning or replacement.
- Base available in #20 Pearl.
- All materials are fire retardant.
- Choose from 13 standard vinyl colours and endless custom selections.

Tech Specs

- Low base height: 53.3 cm lowered, 82.6 cm raised. Base: 56.2 cm by 61.6 cm. Depth upright: 110.5 cm. Reclined: 151.1 cm. Seat: 54.6 cm wide by 38.1 cm deep. Armrest distance: 54.9 cm. Shipping weight: 140.6 kg. Standard electrical: 120 volts, 60 Hz, 3.0 amperes. Optional: 230 volts, 50 Hz, 3.0 amperes.

RA070652019/GT

 HAAG-STREIT
INTERNATIONAL

 **DEVICE
TECHNOLOGIES**

📞 1300 DEVICE (338 423) ✉ customers@device.com.au 🌐 www.device.com.au