

## Assistant's Page: Impression Techniques for MAPA™



Tip of tongue touching soft palate. Practice biting on second molars (pretend if not present).

Goal is to get condyles in fossae to rotate in correct orthopedic position.



With tip of tongue on soft palate, bite as above on cotton rolls, 4 folded gauzes (2x2's) or most commonly on 2 folded gauzes placed on molars (keep off 2nd bicuspids). Goal is to record bite registration just inside freeway space.

# Assistant's Page: Impression Techniques for MAPA™



Using a stable bite registration material,  
record from 2nd bicuspid around to  
opposite 2nd bicuspid.

Contact Pittman Dental Laboratory for more information on  
fabrication of the MAPA

1-800-235-4720

[support@pittmandental.com](mailto:support@pittmandental.com)

[www.pittmandental.com](http://www.pittmandental.com)

Pittman Dental Laboratory is a state of the art  
full service dental laboratory.

233 Centennial Circle

Gainesville, GA 30504

Danny K Crout, DMD, MS [MAPAbiteguard.com](http://MAPAbiteguard.com)