



cm

inch

Place the base of the hearing aid at zero and measure to the top of the curve. Refer to the Fitting Guide on reverse to determine the correct model of Ear Gear for your hearing aid.

FITTING GUIDE

EAR GEAR MICRO

up to 1 in / up to 2.5 cm

EAR GEAR MINI

1 to 1.25 in / 2.5 to 3.2 cm

EAR GEAR ORIGINAL

1.25 - 2 in / 3.2 - 5 cm

EAR GEAR FM

2-3 in / 5 - 7.6 cm

NO MEASURING REQUIRED

COCHLEAR, RONDO & BAHA

Generic model for most sizes

ITE

Clasp and hook security system



250.766.1838

1.888.766.1838

info@gearforears.com

www.gearforears.com