



FULLY CHEMICAL SYNTHESIZED BLOCKING REAGENT

The Blockmaster™ series consists of several synthetic water soluble polymer blocking reagents that solve the safety, quality and purity issues often observed when using “standard” protein based blocking solutions.

Key Properties

- Fully-synthesized polymers
- Water soluble polymer
- Animal protein-free, virus-free
- High lot-to-lot consistency
- Low non-specific binding (NSB) of proteins

Blockmaster™ Series

Grade	Mode of Attachment	Application	Target Materials
CE510, CE210	Primary • Covalent coupling	• Latex agglutination • Chemiluminescent immunoassay • Immunochromatography • Improvement of particle dispersion	• Latex beads • Magnetic beads • Colloidal gold particles
	Secondary • Physical adsorption	• Dispersion stabilizer of colloidal gold particles and magnetic beads	
DB1130	• Physical adsorption	• Latex agglutination • Chemiluminescent immunoassay	• Latex beads • Magnetic beads
PA1080	• Physical adsorption	Anti-adsorption of protein and cell • ELISA Plate • Immunochromatography • Microfluidics chip • 3D cell culture	• PSt, PP, COP, Glass, PDMS, PVDF, Nitro cellulose, etc.

Grade	Package Volume	Solids Content
CE510 CE210	100 mL, 1000 mL	2%
DB1130	100 mL, 1000 mL	10%
PA1080	100 mL, 1000 mL	1%

Contact Information

JSR Corporation

TEL: 81-3-6218-3557

E-mail: dx@jls.jsr.co.jp

URL: <http://www.jsrlifesciences.com>

Medical & Biological Laboratories Co., Ltd. (Japan)

TEL: 81-81-52-238-1904

E-mail: support@mbl.co.jp

URL: <http://ruo.mbl.co.jp/>

JSR Life Sciences, LLC (North America)

TEL: 1-408-543-8800

E-mail: diagnostics.us@jsrlifesciences.com

URL: <http://www.jsrlifesciences.com>

JSR Micro NV (EU)

TEL: 32-16-832-832

E-mail: diagnostics.eu@jsrlifesciences.com

URL: <http://www.jsrlifesciences.com>

MBL Beijing Biotech Co., Ltd (China)

URL: <http://www.bio-med.com.cn>

Beijing

TEL: 86-10-8070-7015

E-mail: yingqi_sun@mbl-china.cn

Shanghai

TEL: 86-21-6468-1752

E-mail: yongjie_yuan@mbl-china.cn,

zhimin_hu@mbl-china.cn

Guangzhou

TEL: 86-20-2212-2010

E-mail: shumei_zhong@mbl-china.cn

Storage

Please store at 2–8° C. Do not freeze the bottle.

Important Notice

- These products are research use only and not intended for therapeutic or in vitro diagnostics use.
- The specifications of these products may be changed without notice.
- JSR Corporation does not guarantee that these products will be continuously available.

JSR Life Sciences Corporation makes no warranties as to this product including, but not limited to, implied warranties of merchantability or fitness to a particular purpose. Blockmaster™ is a global trademark of JSR Corporation. All rights reserved.



JSR Life Sciences

NORTH AMERICA

JSR Life Sciences

408-543-8800

diagnostics.us@jsrlifesciences.com

ASIA

JSR Life Sciences

+81-3-6218-3557

info@jls.jsr.co.jp

EUROPE

JSR Life Sciences

+32-16-832-832

diagnostics.eu@jsrlifesciences.com