

## Easy to use anywhere.

Fill the control unit with ice and water, attach your wrap, customize treatment settings, and you are ready for recovery. Use the optional battery pack for treatment on the go.



FILL ICE AND WATER



SELECT PROGRAM (OPTIONAL)



ADJUST SETTINGS



START TREATMENT

# FINISH STRONG AFTER INJURY OR SURGERY



To rent Game Ready® hassle-free for your recovery, please talk to your physician or visit [gameready.com](http://gameready.com).

AVANOS | GAME READY®

1800 Sutter Street, Suite 500, Concord, CA 94520

888.426.3732 | +1.510.868.2100 | [gameready.com](http://gameready.com)

GAME READY®

## Get next-level cold and compression therapy.

The Game Ready® System combines a powerful control unit with a complete range of anatomic, dual-action wraps to bring proven cold and intermittent compression therapies together with all-new adjustability, precision, and convenience.

*"Awesome for reducing pain and allowing me to rehab faster and more comfortably."\**

P. Klawsy, Patient



\*Based on individual experience, not a manufacturer's claim, results may vary. \*\*Thermal images and skin temperatures taken immediately after 30 minutes of treatment, all devices applied per manufacturer's instructions under identical conditions.

1. Su EP, Perna M, Boettner F, et al. A prospective, multi-center randomized trial to evaluate the efficacy of a cryopneumatic device on total knee arthroplasty recovery. J Bone & Joint Surg. 2012;SuppA(94-B):153-6. 2. Waterman W, Walker JJ, Swains C, et al. The efficacy of combined cryotherapy and compression compared with cryotherapy alone following anterior cruciate ligament reconstruction. J Knee Surg. 2012;25(02):155-160. 3. Murgier J and Cassard X. Cryotherapy with dynamic intermittent compression for analgesia after anterior cruciate ligament reconstruction. Preliminary study. Orthop Traumatol Surg Res. 2014;100(3):309-312. 4. As with all RICE applications. 5. Leegwater NC, Jore HW, Brohet R, Nolte P. Cryocompression therapy after elective arthroplasty of the hip. Hip Int. 2012;22(05):527-533.

There are inherent risks in all medical devices. Please refer to the product labeling for Indications, Cautions, Warnings and Contraindications. Refer to [www.gameready.com](http://www.gameready.com) for product safety technical bulletins.

GAME READY®

ACTIVATE NEXT-LEVEL INJURY & POST-OP RECOVERY



**Cold Therapy.** Circulating ice water continuously and consistently cools your injury site to safely and comfortably reduce pain and swelling.



**Dual-Action Wraps.** Anatomically designed, circumferential wraps surround your injury site, covering more surface area for deeper cooling.



**Intermittent Compression.** Gently squeezes and releases to remove swelling, stimulate the flow of blood, and conform the wrap to your body.



**Simple Controls.** The system makes it easy to use pre-set programs and quickly customize temperature, compression, and time settings.

# FINISH STRONG AND GET BACK TO BEING YOU

After an orthopedic injury or surgery, you play a critical role in the recovery process. The Game Ready® cold and compression system allows you to reduce pain and swelling while reducing the use of opioid medications. All so you can finish strong, maximize your rehab, and get back to the things that matter as quickly as possible.

## Faster, deeper, longer-lasting cooling for proven therapeutic benefits:

- ✓ Reduces your opioid consumption and may allow for earlier discontinuation of opioid medications in certain patients.<sup>1,2,3</sup>
- ✓ Decreases pain and swelling.<sup>1,4</sup>
- ✓ Can increase post-operative range of motion.<sup>3</sup>
- ✓ May increase your satisfaction with the recovery process.<sup>1,5</sup>

*“Game Ready provides a significant decrease in the need for pain medication and at least a 20% faster recovery.”\**

Peter Millett, MD, MSC | The Steadman Clinic

You play a critical role in your recovery after injury or surgery. Game Ready is here to help.



## Innovative wrap design for better comfort and results.

Anatomically engineered for all major body parts, Game Ready's patented wraps ensure better fit, greater comfort, and circumferential coverage for more effective therapy.



Trusted by top physicians, professional sports teams, elite athletes, and special military forces around the world, the Game Ready System takes traditional RICE methods (Rest-Ice-Compression-Elevation) to the next level.

## THERMAL IMAGERY SHOWS THAT GAME READY OUT-COOLS THE COMPETITION\*\*

