

FROM BULK-FED TUBES OR VIALS TO PERFECTLY LABELED—IN MINUTES

FLEX TUBE PRO

The more automated a high-volume lab becomes, the more productive it will be. And automating the positioning and labeling steps for tubes and vials will increase throughput. It might even improve staff morale.

Our new Flex Tube Pro takes the reliability of the Tube Pro PR100i printer-applier and adds a front-end bulk feeder to improve labeling throughput. Up to 2000 tubes or vials can be placed in the gravity-fed hopper at a time. The system will automatically orient and position them, print a flawless barcode label, and apply it in exactly the same location on the tube, time after time.

FIVE MAJOR FEATURES DISTINGUISH THIS SYSTEM FROM ANYTHING YOU'VE SEEN BEFORE:

- It's specifically designed to wrap labels around cylindrical objects, such as tubes and vials;
- Label placement is consistent and accurate;
- High resolution barcode printing provides outstanding read rates;
- Quick changeover makes it easy to accommodate tubes of multiple lengths and diameters;
- And an array of interface options makes it easy to integrate with the LIS and other automation in the lab.

Tubes and vials can be plastic or tempered glass, either capped or uncapped. Vials of any outside diameter from 6 mm to 15 mm can be handled, and they can measure from 12 mm to 90 mm in length. Placement accuracy is ± 0.2 mm; barcodes can be linear or 2D and are printed in high resolution with a 600 dpi (dots per inch) printer.



YOUR NEEDS, OUR EXPERTISE.

Come to us with your laboratory automation and sample identification requirements. We will create a solution that not only meets your needs but will exceed your expectations.

LOOK TO US FOR ALL YOUR CRITICAL LABELING SOLUTIONS.

We work closely with you to provide the right labels to meet all your requirements. And not just the right labels, but everything for your labeling needs. From labels to printers, and scanners to label replication systems, Computype offers this and everything in between.

WE'D LOVE TO HEAR FROM YOU.

Computype
2285 West County Road C
St. Paul, MN 55113-2567
Phone: 800-328-0852 | 651-633-0633
Fax: 651-633-5580
Email: sales@computype.com
Twitter: [@computype](https://twitter.com/computype)

www.computype.com

Computype
Kingswood Business Park, Connaught Road
Kingston Upon Hull
HU7 3AP
United Kingdom
Phone: [44] 1482-835366
Fax: [44] 1482-831152 421-606
Email: contactus@computype.com

www.computype.com/eng-gb

LS0924.r1.1013