

# **Cryoport Express® Shipper Accessories Specifications**

Cryoport offers a variety of accessories available to be included in your shipment, including cryoboxes, cryovials, canes and Tyvek bags. If accessories are needed for your shipment, you can select them at the time of your order or create standing orders with pre-arranged accessory configurations through the Cryoportal™, our logistics management platform.

## **Cryoboxes:**

A cryobox is a storage box for frozen materials. They are used for storage and shipment of cryovials. The boxes can be shipped flat in the dewar box and assembled on site.

Material:	Cardboard
Dimensions:	2" x 2" x 5" (50.8 x 50.8 x 127 mm)
Storage Capacity:	Thirty 2.0 ml cryovials

### **Cryovials:**

Cryovials are vials used for storing biological materials at frozen temperatures.

Material:	Polypropylene
Dimensions:	2.0 ml with graduations

#### Canes:

Canes are used to hold cryovials in liquid nitrogen containers.

Material:	Aluminum
Dimensions:	.5" diameter x 11.25" tall (12.7 x 285.75 mm)
Storage Capacity:	Five 2.0 ml cryovials per cane

#### Tyvek Bags:

Tyvek bags are used for infectious substances and pass IATA testing for Class 6.2 substances.

Material:	Tyvek – 1073B (medical grade)
Dimensions:	5.5" wide x 18" tall (folded 11" tall) (139.7 x 457.2 mm)
Storage Capacity:	15 canes or 2 cryovial boxes