

FOR USE WITH AIRX FRESHENAIRES & SCENTINELS

EVEN WITH COMPETITIVE CABINETS



RX 19E MICROCELL ODOR COUNTERACTANTS

RX 19E microcells are the nearest thing to a liquid in solid form. They hold practically 100% of their calculable volume of AIRICIDE® Odor Counteractant, yet they won't spill or drip like liquids. Nor will they skin or crust over as do gels. The space age microcells hold the AIRICIDE® till exposed to the air, then gradually, but steadily, release it by capillary action to change the size and shape of odorous molecules so it is

recognized as something different and more pleasant that it would have been.

Odors float in the air in the form of molecules that when sniffed in are identified by minute hair-like antenna...The microcell blocks are small and compact (2½x4x½) allowing use of multiple blocks in a single circulator or wall cabinet and adapters available allow use of 19E with just about any competitive cabinet.



Block is easily removed from poly bag without even touching block.



Use RX 19Es with Airx Freshenaire to control intense foul odors.




Use adapters in existing and competitive cabinets.



Scentinel Jr. controls release thru adjustable fins. No battery or fans.

AIRICIDE®, RX 19E MICROCELLS & CABINETS



1. AIRICIDE®

Not a mask, not a desensitizer, AIRICIDE® is formulated to eliminate foul odors by changing the geometrics of odorous molecules in the air, so when they reach cilia receptors in the nose they are not the bad odors they would have been.

2. MICROCELLS

A space age, rigid polymer foam that holds nominally 100% of its dimensional (2½"x4"x½") capacity. Then gradually and constantly releases its vapors into the air to neutralize foul odors from cancer, burns, human incontenance, animals, decay, bacterial growth, fire, smoke, liquor, beer, mildew, etc.

SPECIFICATIONS

3FECIFICATIONS	
TYPEc	ompressed microcells
ODOR COUNTERACTANT	pure AIRICIDE®
RELEASE OF COUNTERACTANT	approx. 60-90 days
OLFACTORY DESENSITIZERS	
PETROLEUM DISTILLATES	none
TOXIC	no
EFFECT ON METALS	non-corrosive
COLOR	white, non-staining
FREEZE-THAW STABILITY	no degradation
HEAT STABILITY	no degradation
SHELF LIFE	min. 1 year, closed
DRYING CHARACTERISTICSnc	shrinking, no skinning
VERSATILITYneutralizes odors oxidation re	from bacterial growth, eduction, chemicals, etc.

3. THE AIRX FRESHENAIRE

A portable case (13"x16"x4") with handle. Built-in extension cord and 110v electric fan to circulate *AIRICIDE**, and cover to conserve microcells when not being used.

5. THE SCENTINEL WALL CABINET

when installing a new block.

Brushed metallic cabinet with computerized mechanics operated by a single 'D' battery. Light on front signals low battery. Can be programmed to turn off when room is not occupied. Use with 1,2 or 3 microcells, depending on need.

These foam pieces allow conversion of any deodor-

izer cabinet to the use of RX 19E Microcells.

Competitive cabinets are more effective with


AIRICIDE®, and adapters can hold two microcells.

There's no need to discard one that still has life

6. THE SCENTINEL JR.

The simple solution to many odor problems. Holds a single RX 19 microcell. No battery to replace or be stolen. No fan – Fins on side of cabinet allow adjustment for individual needs.

DISTRIBUTED BY:



3 Chattanooga Irvine, CA 92620-2514

Tel 714-832-8441

chemex

info@chemexindustries.com