

# LabMinder

## Laboratory Monitoring System

**AND**  
Discover Precision

### 24/7/365 Unattended Operation with Automated Notification and Response


Today's lab manager is under constant pressure to increase productivity, but at the same time maintain or even reduce costs. One of the quickest ways to increase output is by simply adding operating hours. But without any other changes, the benefit would quickly be offset by the cost of the additional manpower to support those hours.

LabMinder is a stand-alone lab-wide system that supports unattended operation by providing automated

notification and response capabilities when no operator is present. Unmanned test systems are continuously monitored, and designated personnel are notified by text message or email if an important event is detected. Responses can be sent back to LabMinder, triggering predefined actions in case of error.

Notification can either be sent to an entire list simultaneously, or to the first person on the list, and then escalated through the list until someone responds.

*LabMinder constantly monitors and interacts with the various systems within the lab, and notifies designated personnel via text or email if an important event is detected.*




### Benefits

- Reduced overhead requirements with unattended operation capabilities
- Increased efficiency through automated notification and response
- Maximized testing time through prompt fault detection and correction

### Features


- Centralized monitoring of all test systems
- Choice of SMS, email or push notifications
- Health metric customizable by lab personnel; allows monitoring of any condition in a running system
- Compatible with most existing mobile apps, including email, text and web messaging
- Individual system override capability

## LabMinder SMS Notifications



LabMinder typically sits inside the lab rather than on a company network, and is connected directly to an SMS modem.

## LabMinder Email Notifications



LabMinder sends email via an existing SMTP server. The response is displayed via a temporary web page on the LabMinder display.

### Americas

A&D Technology, Inc.  
Ann Arbor, MI USA  
Ph: +1 (734)973 1111  
www.AandDTech.com

### Europe

A&D Europe, GmbH  
Griesheim, Germany  
PH: +49 (0)6155-605 250  
www.AandDEurope.com

A&D Europe - UK Branch  
Abingdon, Oxon, UK  
PH: +44 (0)1235-550 420  
www.AandDEurope.com

### Asia

A&D Technology Trading Co.  
Shanghai, China  
PH: +86 (0)21-3393 2340  
www.AandDTech.com

A&D Company, Ltd.  
Tokyo, Japan  
PH: +81 (0)3-5391 2753  
www.AandD.co.jp

### Australia

A&D Australasia Pty Ltd.  
Kensington, VIC, Australia  
PH: +61 (0)3-9372 1522  
www.andaustralasia.com.au