

# Different Sizes. All Smart. RFID-Enabled Cabinets.

Hospital pharmacies are under increasing pressure to conserve costs while maintaining adequate and safe inventory. The IntelliGuard® Smart Cabinet Inventory System™ uses RFID technology to provide real-time visibility and reduce hospital inventory costs for optimal management of high-value, critical dose pharmaceuticals.

## Smart Cabinet Inventory System<sup>™</sup> (SCIS)

IntelliGuard® Smart Cabinets<sup>™</sup> monitor critical drug levels and temperature storage conditions. Inventory levels are automatically updated as items are removed or replaced. That data is relayed to the IntelliGuard® Client Intelligence Portal for instant access to actionable inventory information in real-time.

The IntelliGuard<sup>®</sup> Smart Cabinets<sup>™</sup> are available in multiple sizes — from small under counter to large capacity models — to fit your facility.


#### **RFID Technology Puts You In Control**

- Increased security with access control and real-time visibility to user transactions
- Confidently reduce inventory and pre-paid inventory costs
- · Improve patient safety by reducing stock-outs
- Manage inventory expiry details, storage temperatures and excursions at the item level
- Enable truly effective inventory control from the web-based IntelliGuard®
  Client Intelligence Portal

Actionable alerts and real-time visibility in the palm of your hand.



## **HOW INTELLIGUARD® SCIS<sup>™</sup> WORKS**



- Open the Door: Secure Remove or Restock Badge Log-in Optional Items
- Close the Door
- View Real-Time Inventory via Web Portal

## Track Critical Medication Inventories Automatically.

## The IntelliGuard<sup>®</sup> Smart Cabinets™

IntelliGuard provides one solution to store and dispense Cold (2 - 8°C) as well as Controlled Room (20 - 25°C) temperature medication.

Using the latest RFID technology, IntelliGuard provides real-time visibility of your pharmaceutical stock. Quickly find and remove expired or recalled medication, and identify shortages to avoid stock-outs that put patients at risk.

### **Intelligent Inventory Management**

The IntelliGuard<sup>®</sup> Smart Cabinets<sup>™</sup> keep you connected to your inventory through actionable alerts and reports.

l	Guard			
intein	Guard	Dashboard	Dashboard	
	Deshboard	Cabinet Activity	> Welcome System Administrator,	
any in	Welcome System Administrator, Below is a summary drug must neare dom and activity. Use the menu to access detailed information about status of the subole and more threation.	Alerts	Below is a summary of your most recent alerts and activity. Use the menu to access detailed information about the	
teriteda a	Peter relation to determine address	Reports	status of the cabinet and your inventory.	
olart Vachar 3		Inventory Levels	> Pacility cedition V Cabinet cedition V	
	Kets - N Facility - Al Cabrels	inet Tracker	> Tenty COLLY	
	Protect Receille (4)  Steenbary MALLevels (5)		Alerts – All Facilities – All Cabinets	
	Storf-Genet Products (2)		Product Recalls (4)	
	Recert Activity - Al Facilies - Al Carries			
	There is an admity for today		🤊 Inventory PAR Levels (0) 🔤 👘	
	O There is no activity for yearentary			
	Approved Inventory - All facilies		5 Short-dated Products (2)	
	Product forcats (4)			
	Product Terrate (4)		Recent Activity – All Facilities – All Cabinets	
	Standary Millionnia (S)			
	Shan-baket Products (2)			



## Put IntelliGuard® to work for you.

Contact us today for a live demo of the IntelliGuard® SCIS™. Call 760-448-9500 or visit http://IG.Solutions.



- Wi-Fi or 3G/4G
- Secure access control panel
- Real-time temperature monitoring
- Adjustable and fixed shelves accommodate from 50 to over 500 items
- **Critical Inventory Alerts** (par level, expiry)

© 2018 MEPS Real-Time, Inc. All rights reserved. IntelliGuard® is a dba of MEPS Real-Time, Inc. MEPS, MEPS Real-Time and IntelliGuard are registered trademarks and Intelligent Inventory Solutions, IG Solutions and the IG logo are trademarks of MEPS Real-Time, Inc., Carlsbad, CA, in the United States and/or other countries. 900-0101. Rev D