

EAU CONTOU

808 Clay

Warm

554 3.0 Areola 3.0

Pigments



p823155813

822 Arab Brown

100 Just Black

The NEW Nouveau Contour pigments:

MEAU CO

more colors, more choice, for an even more natural result.

Organic pigment line 🕟

Very high concentration of pigment, bright colors, limited follow-up treatments

The Organic pigments of Nouveau Contour are organic and free of iron oxides: your guarantee that they are void of nickel impurities. The product composition contains less alcohol: the same pigment volume contains more pigment. The result? A higher concentration of pigment will remain in the skin, resulting in bright colors. An additional advantage is that the low concentration of alcohol will dry out the pigment less quickly.

The formulation is developed with unique bonding materials, allowing for optimal skin implementation and pigment retention. The Organic pigments have excellent flow properties which enables effortless pigment deposit. During the procedure, the pigment is more easily inserted into the skin. The colors remain clear and fresh. After healing, the results are very natural and follow-up treatments are limited.

The Organic pigment range includes 38 colors, 9 eyebrow colors, 3 eye colors, 12 lip colors, 5 areola colors, 4 camouflage colors and 5 areola colors.



• 100% organic pigments

CO

202

Ginger Bread

Warm

- Free from iron oxides
- Less alcohol
- Higher concentration of pigment
- Bright colors
- Very natural results







Eyeliner pigments





212 - Deep Brown





104 - Admiral



214 - Black Cool Brown







216 - Chocolate

Camouflage pigments





457 - Camo D.O



453 - Camo B.O



455 - Camo C.O

Areola pigments



556 - Areola 4.0





558 - Areola 5.0



554 - Areola 3.0

Important note: the pigment colors in this brochure are for illustrative purposes only and may differ from reality. For a true-to-color representation of each pigment color, please refer to the CIS guide or the actual pigment.

Organic pigment line

Correction pigments





Note: These results are achieved with Nouveau Contour pigments and shown for illustrative purposes. These are not a result of the new Nouveau Contour pigments. Do you want to share your work? Send us your pics: **social@nouveaucontour.com**

Photocredits: Magdalena Mazurek, Hanna Idczak Nouveau Contour Poland Partner



Fusion pigment line 🛞

Bright, intense colors and excellent pigment retention

The Fusion line combines organic and inorganic ingredients in an innovative pigment formula.

The advantage of the inorganic component is that the pigments remain very stable in the skin. The safe iron-oxides meet all stringent international regulations regarding nickel impurities. These carefully selected iron-oxides will give you the extra boost of nearly immediate real life color in the skin, intense color saturation and optimal color retention, therefore nearly eliminating touch-up procedures.

The organic component -free from iron oxides- ensures the clear color and good pigment retention. After healing, the colors are bright, extremely soft and very natural.

This pigment line is indispensable for eyebrow treatments. Nouveau Contour Fusion pigments are ideal for the ombre technique with amazing results. Hair strokes can also be set very accurately and perfectly.

The Fusion line includes 12 eyebrow colors, 10 lip colors and 4 correction colors.

- Innovative pigment formula
- Exceptional healed results
- Very natural results with minimal touch-ups
- Excellent pigment retention
- Brilliant & intense color tones
- Eyebrow colors with balanced formula to offset grey or red healing in the skin
- Lip colors do not change: they only get slightly lighter



Eyebrow pigments



888 - Soft Mocca



Color Identification System

Nouveau Contour's unique Color Identification System simplifies the complexity of color theory, making it easy to choose the correct pigment for each procedure and every skin type. All of the relevant information is available for each pigment color, providing technicians with added confidence in the color outcome of the skin.

The Color Identification System helps you to:

- Choose the correct pigment for each procedure & every skin type
- Immediately identify which pigment to add to make a color darker or lighter
- Recognize if the undertone is warm, cool or neutral
- Obtain a clear description of the pigment

- UNIQUE!
- Simplifies color selection
- Choose the correct pigment for each procedure & every skin type





MASTERS IN MICROPIGMENTATION

Pigment display

Show pigments in style on this universal pigment display with LED. The display can show up to 18 pigment bottles. Next to a general top card, there are 4 topcards of choice: general, lips, eyes and brows.

ALARY RATE RATE RATE RATE

Contact us! 111 N. Orange Ave, Suite #775 | Orlando, Florida 32801 P. Ph: (321) 332-6975 Toll Free: (877) 628.7201 E. info@nouveaucontourusa.com | www.nouveaucontourusa.com