

Light Up Your Research With Fluorescent Proteins

Explore the entire fluorescent protein plasmid collection at addgene.org/fluorescent-proteins



A better way to share plasmids

Join the Addgene Community:


[@Addgene](#) [Addgene](#) [company/Addgene](#)

www.addgene.org

help@addgene.org


Key:
Fluorophore
excitation / emission
wavelength (nm)
*LSS- Long Stokes Shift

Excitation and Emission Wavelengths of Sample Fluorescent Proteins



Subcellular Localization

These are representative plasmids. For a more exhaustive list, please visit addgene.org/fluorescent-proteins/localization



Special thanks to the Michael Davidson Lab for providing fluorescent images. Browse 3000+ plasmids from the Davidson Collection at: addgene.org/fluorescent-proteins/davidson