

THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY™



 **thermoflect**® Heat Reflective Technology™

Today's fast-paced surgical environment needs products that work, not products that make work.

Thermoflect heat reflective technology products are simple, safe and effective. They begin working from the moment they are applied to bank heat and help prevent redistribution temperature drop [RTD]. They use science instead of electricity and are the most efficient and cost-effective prewarming products available.

We invite you to try Thermoflect heat reflective technology products in your facility. Learn more at www.Thermoflect.com.



**NOT FOR
USE IN MRI**

 **encompass**®

Better care starts
with safety and comfort.

THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY™



Thermoflect® Patient Care Products

CAT NO	DESCRIPTION	SIZE	PACKAGING	QTY/CASE
5100-200	Blanket, Small	36" x 40"	5 EA/BG, 12 BG/CS	60 EA
5100-400	Blanket, Medium	48" x 48"	5 EA/BG, 10 BG/CS	50 EA
5100-700	Blanket, Large	48" x 84"	5 EA/BG, 5 BG/CS	25 EA
5120-350	Leg Coverings	34" Long	5 PR/BG, 10 BG/CS	50 PR
5200-400	Procedure Pak 400 Medium Blanket, Leg Coverings, Adult Cap	Universal	Individually Packaged	20 EA
5200-700	Procedure Pak 700 Large Blanket, Adult Cap	Universal	Individually Packaged	24 EA
5105-100	Pulse Oximeter Mitt	Universal	5 EA/BG, 10 BG/CS	50 EA
5300-100	TherMuff™ Bag, Sterile	8.5" x 9"	Individually Packaged	50 EA

Thermoflect® Staff Apparel

CAT NO	DESCRIPTION	SIZE	PACKAGING	QTY/CASE
5140-100	Jacket, Silver Side Out	X-Small	5 EA/BG, 5 BG/CS	25 EA
5140-150	Jacket, Blue Side Out	X-Small	5 EA/BG, 5 BG/CS	25 EA
5140-200	Jacket, Silver Side Out	Small	5 EA/BG, 5 BG/CS	25 EA
5140-250	Jacket, Blue Side Out	Small	5 EA/BG, 5 BG/CS	25 EA
5140-300	Jacket, Silver Side Out	Medium	5 EA/BG, 5 BG/CS	25 EA
5140-350	Jacket, Blue Side Out	Medium	5 EA/BG, 5 BG/CS	25 EA
5140-400	Jacket, Silver Side Out	Large	5 EA/BG, 5 BG/CS	25 EA
5140-450	Jacket, Blue Side Out	Large	5 EA/BG, 5 BG/CS	25 EA
5140-500	Jacket, Silver Side Out	X-Large	5 EA/BG, 5 BG/CS	25 EA
5140-550	Jacket, Blue Side Out	X-Large	5 EA/BG, 5 BG/CS	25 EA
5130-200	Vest, Reversible	Small	5 EA/BG, 5 BG/CS	25 EA
5130-300	Vest, Reversible	Medium	5 EA/BG, 5 BG/CS	25 EA
5130-400	Vest, Reversible	Large	5 EA/BG, 5 BG/CS	25 EA
5130-500	Vest, Reversible	X-Large	5 EA/BG, 5 BG/CS	25 EA



Thermoflect® Trauma Products

CAT NO	DESCRIPTION	SIZE	PACKAGING	QTY/CASE
5170-200	Transport Cocoon™	Adult	Individually Vacuum Packaged	10 EA
5170-150	Transport Cocoon™	Pediatric	Individually Vacuum Packaged	10 EA
5100-250	Papoose	5 to 10 lbs	5 EA/BG, 10 BG/CS	50 EA



TECHstyles®

Superior Nonwoven Products for Infection Control



All items featured on this page are free of natural rubber latex

800-826-4490 | www.thermoflect.com

© 2013 Encompass Group, LLC
All Rights Reserved
Lit. No. AC10145 rev. 07/13



5145-500 Standard Gown



5145-500IV IV Access Gown



5145-500M Total Access Gown
(More information on back panel)

Thermoflect® Adult Patient Apparel and Prewarming Admission Kits

CAT NO	DESCRIPTION	SIZE	PACKAGING	QTY/CASE
5110-100	Bouffant Cap	Adult	5 EA/BG, 20 BG/CS	100 EA
5145-500	Standard Gown	Universal	5 EA/BG, 5 BG/CS	25 EA
5145-500IV	IV Access Gown	X-Large	5 EA/BG, 5 BG/CS	25 EA
5190-600	IV Access Gown	Small	5 EA/BG, 5 BG/CS	25 EA
5145-500M	Total Access Gown	X-Large	5 EA/BG, 5 BG/CS	25 EA
5115-200	Footies	Universal	5 PR/BG, 5 BG/CS	25 PR
5000-800	Admission Kit, IV Access Gown, Adult Cap, Patient Safety Footwear, Belongings Bag	X-Large	Individually Packaged	25 EA
5000-810	Admission Kit, IV Access Gown, Adult Cap, Patient Safety Footwear, Belongings Bags (2)	X-Large	Individually Packaged	25 EA
5000-820	Admission Kit, IV Access Gown, Adult Cap, Patient Safety Footwear	X-Large	Individually Packaged	25 EA
5000-830	Admission Kit, IV Access Gown, Patient Safety Footwear	X-Large	Individually Packaged	25 EA
5000-840	Admission Kit, IV Access Gown, Adult Cap	X-Large	Individually Packaged	25 EA
5000-850	Admission Kit, IV Access Gown, Patient Safety Footwear, Belongings Bag	X-Large	Individually Packaged	25 EA



Thermoflect® Pediatric Patient Apparel

CAT NO	DESCRIPTION	SIZE	PACKAGING	QTY/CASE
5110-250	Cap With Stockinette Comfort Band	Neonatal	5 EA/BG, 20 BG/CS	100 EA
5110-200	Cap With Stockinette Comfort Band	Infant	5 EA/BG, 20 BG/CS	100 EA
5110-150	Cap, Bouffant	Youth	5 EA/BG, 10 BG/CS	50 EA
5180-200	Jacket (2-4)	X-Small	5 EA/BG, 10 BG/CS	50 EA
5180-300	Jacket (6-8)	Small	5 EA/BG, 10 BG/CS	50 EA
5180-400	Jacket (10-12)	Medium	5 EA/BG, 10 BG/CS	50 EA
5180-500	Jacket (14-16)	Large	5 EA/BG, 10 BG/CS	50 EA
5185-200	Pants (2-4)	X-Small	5 EA/BG, 10 BG/CS	50 EA
5185-300	Pants (6-8)	Small	5 EA/BG, 10 BG/CS	50 EA
5185-400	Pants (10-12)	Medium	5 EA/BG, 10 BG/CS	50 EA
5185-500	Pants (14-16)	Large	5 EA/BG, 10 BG/CS	50 EA
5190-200	Pediatric Gown (2-4)	Small	5 EA/BG, 5 BG/CS	25 EA
5190-300	Pediatric Gown (6-8)	Medium	5 EA/BG, 5 BG/CS	25 EA
5190-500	Pediatric Gown (10-12)	Large	5 EA/BG, 5 BG/CS	25 EA



THERMOFLECT® HEAT REFLECTIVE TECHNOLOGY™ TOTAL ACCESS PATIENT WARMING GOWN



ThermoFlect Heat Reflective Technology Total Access Gown™

The ThermoFlect Total Access Gown is designed to provide maximum access to the surgical patient without compromising the benefits of patient warming. The unique 360 degree release strip provides quick separation for access to the surgical site and simple repositioning for transport to PACU.

This multi-function gown can be easily configured in the four positions below as well as a multitude of others to address the specific needs of procedural positioning.

- **Full Length** for prewarming and recovery
- **Open Arm** for chest or breast surgery
- **Upper Body Blanket** for abdominal or lower body procedures
- **Lower Body Blanket** for abdominal or upper body procedures

Evidence-Based Practices for ThermoFlect Heat Reflective Technology

The Effect Of Preoperative Reflective Hats And Jackets, And Intraoperative Reflective Blankets On Perioperative Temperature, The Internet Journal of Anesthesiology. 2003. Volume 6, Yiwei Sheng, Frank Zavisca, Elizabeth Schonlau, Renee Desmarattes, Edwin Herron, Randall Cork, Department of Anesthesiology, LSUHSC Shreveport

Assessing the Use of a Reflective Blanket in Managing Post Operative Hypothermia, Rush University Medical Center and Rush University College of Nursing, Chicago IL, Kimberly Humbarger BSN RN, Isabel Pecson RN, Ruth Kleinpell PhD RN FAAN FCCM

The Usefulness of Preoperative Warming, B. van Beek, CRNA, Anaesthesia Department - Queen Beatrix Hospital, Winterswijk, The Netherlands, TSCI Abstract, October 2010

Thermal Efficiency of Pre-warmed Cotton, Reflective, and Forced Warm Air Inflatable Blankets in Trauma Patients, Emergency Nurses Association, 2000 Conference, Research Paper Presentation

Maintenance of Normothermia During Relatively Minor Orthopedic Surgical Procedures On the Extremities, Leo Leliveld, CRNA, Best Practices Lit. No. AC11187, January 2011

A Comparative Study of The Thermal Efficiency Of Thermo-Lite Clothing Versus Pre-Warmed Cotton Blankets In Preventing Intraoperative Hypothermia, By Joanne A. Kelley R.N., B.S.

