



SUGGESTED FORMULA

TRICHLORACETIC ACID 35% (W/V) TOPICAL SOLUTION

Version number: 1.0

Volume: 100 ML

Trichloroacetic Acid ACS Reagent Crystals	35 g
Purified Water, USP	QS 100 mL

SUGGESTED COMPOUNDING PROCEDURES

1. Add Trichloroacetic acid to purified water in calibrated, glass bottle equal to 50% of the final volume.
2. Bring formula to final volume with purified water, cap shake to mix until uniform.
3. Suggested quality assessments:
 - a. Clarity
 - b. Color
 - c. Label (s)- auxiliary labels (external use), storage, BUD

Protect from light.

Store in amber glass bottle

Store at Controlled Room Temperature

No claims are made as to the safety or efficacy of this preparation. This formulation is provided solely at the unsolicited request of the pharmacist.

Beyond-Use Dates of preparations are conservative estimates by the formulator using reference books, peer-reviewed literature, intended duration of therapy, formulation from commercially available products, organoleptic observations and current USP guidelines. Compounders may have stability studies performed by a reputable laboratory if they wish to extend the Beyond-Use Date. It is recommended that you follow USP <795> recommendations for potency testing.

Beyond-Use Date is estimated to be 90 days.

Precautions should be taken to prevent cross-contamination and exposure of ingredients to the compounder and contamination of the preparation by the compounder. Wear appropriate protective equipment . Use safety enclosures (hoods) when weighing and mixing.

Phone: 800.370.6231 | Fax: 732.608.5420 | Internet: SpectrumRx.com | Email: sales@spectrumrx.com
© 2015 Spectrum Pharmacy Products. All rights reserved.

